



E-spares cat Now Available Online
at www.invacare.co.uk

Customer Name :	Account No :
Dealer :	Order No :
Address :	Order Qty :
Tel No:	Order Date :
Delivery address :	

Invacare® TDX SP
DEALER
Max User Weight 180Kg (28 Stone)

Standard Features Standard Seating Manual Tilt In Space Manual Recliner Back Multi Adjustable Legrest ACS 2 Control System Suspension Sure Step 73 Ah Batteries Headrest 8km/h 4 Pole Motors Base Price £4,195	
--	--

Please make your selection with a cross in the box ◆ Standard feature ○ Free of Charge
NC = Incompatibility with option XXX99

Customer weight =		Kg	
Code	<input checked="" type="checkbox"/>	Description	Standard Dealer
ARA0000	<input checked="" type="checkbox"/>	TDX Basic Model	£4,195
SEATING TYPE			
ARA1150	<input checked="" type="checkbox"/>	Standard Seat Unit	◆
BACK CUSHION			
ARA1178	<input checked="" type="checkbox"/>	Back Sling tension adjustable	◆
UPHOLSTERY COLOUR			
ARA0268	<input checked="" type="checkbox"/>	Black (Nylon)	◆
SEAT WIDTH 44-48cm & SEAT CUSHION			
ARA0343	<input checked="" type="checkbox"/>	44 - 48cm	◆
ARA1336	<input checked="" type="checkbox"/>	Seat Cushion (Width 43cm 5cm thick)	◆
SEAT DEPTH			
ARA0326		41 cm (NC 0953, 0954, 0947)	○
ARA0327		46 cm	○
ARA0330		51 cm	○
BACKREST HEIGHT (Without Seat Cushion)			
ARA0448		48 cm	○
ARA0454		54 cm	○
ARMREST			
ARA0635		Fixed Armrest 25 - 34cm height adjustable, removable	○
ARA0636		Fixed Armrest 29 - 38cm height adjustable, removable	○
ARMPAD			
ARA0647		Standard Armpad (32,5 cm)	◆
ARA0649		Armpad extended (+ 5cm) (NC 0602, 0603)	○
SEAT ACCESSORIES			
ARA1353	<input checked="" type="checkbox"/>	Contour Headrest	◆

Order No:

Customer Ref:

LPF0099Q
TDX SP S

Subject to Invacares standard terms conditions.
Invacare reserve the right to alter product specification without prior notice.

SEAT ADJUSTMENTS			
SEAT LIFTER/TILT			
ARA1002		Manually settable Tilt (max user weight 180 kgs)	○
ARA1014		Powered Tilt (up to 28°) (max. user weight 160 kg) NC 0556, 0557	£808
ARA1021		Seat Lifter with electrical Tilt (max user weight 160 kg) 30cm rise 28° tilt NC 1639, 0556, 0557, 1016	£1,540
BACKREST RECLINER			
ARA0958		Settable Recliner (-10° / 0° / 7,5° / 15° / 22,5° / 30°)	○
ARA0953		Pneumatic Recliner RH NC (0326)	£155
ARA0954		Pneumatic recliner LH NC (0326)	£155
ARA0947		Electric Recliner NC (0326)	£360
PLEASE CHOOSE EITHER CENTRE MOUNTED OR SIDE MOUNTED LEGRESTS			
CENTRE MOUNTED LEGRESTS			
FOOTREST HANGER TYPE (INCLUDING CALF PAD AND ANGLE ADJ FOOTPLATE) KNEE TO HEEL 310mm to 410mm			
ARA0556		Hanger , centre mounted, angle settable NC (1014, 1021)	○
ARA0557		Hanger, centre mounted, length compensating, electric. driven NC (1014, 1021)	£420
CALF PADS			
ARA0583		Calf pads	○
FOOTPLATE TYPE			
ARA0559		Angle Adjustable Footplate	○
SIDE MOUNTED LEGRESTS			
LEGREST FIXATION			
ARA0590		Position for seat cushion up to 7cm thick	◆
FOOTREST HANGER TYPE; PAIRS (INCLUDING CALF PAD)			
ARA0503		Hanger Vari F	○
ARA0533		Hanger Vari A NC (0512)	£87
ARA0532		Hanger ADM 0° - 80° NC (0512)	£155
ARA0531		Hanger ADE 0° - 80° NC (0512)	£455
FOOTPLATE TYPE; PAIRS			
ARA0512		Fixed Footplate	○
ARA0504		Angle Adjustable Footplate	£51
ARA0513		Angle and depth adjustable Footplate (Adapt)	£82
REMOTE CONTROL			
USER CONTROL + BRACKET			
ARA1639		REM Type B (up to 4 electric actuators and light functionality)	◆
ARA1617		REM 550 with LCD screen	£82
ARA1601		ACS REM24SD	£31
ARA1614		G90A	£31
ARA1615		G91S Includes goose neck bracketry, select option 1207 or 1208 (NC 1707, 1708, 1609)	£825
		Head-Remote, ultrasonic proximity sensor (only available as special) Please contact UK Specials Dept	POA
ARA1207		Standard Bracket RH (NC 1404)	○
ARA1208		Standard Bracket LH (NC 1405)	○
ARA1707		Swing Away + height adjust. RH High End (NC 1404, 1415)	£97
ARA1708		Swing Away + height adjust. LH High End (NC 1405, 1339)	£97
ATTENDANT REMOTE + BRACKET			
ARA1644		REM Type B (up to 4 electric actuators and light functionality)	£12
ARA1618		REM 550 with LCD screen	£93
ARA1602		ACS REM24SD	£42
ARA1616		G90A	£42
ARA1207		Standard Bracket RH	○
ARA1208		Standard Bracket LH	○
TRAY & CENTRE CONTROLS			
ARA1608		A C S Centre Control (REM 24SD)	£52
ARA1607		A C S Tray Control c/w Acrylic Tray & Actuator Functionality - DX REM 41	£147
ARA1207		Standard Bracket RH	○
ARA1208		Standard Bracket LH	○

Order No:

Customer Ref:

LPF0099Q
TDX SP SSubject to Invacare standard terms conditions.
Invacare reserve the right to alter product specification without prior notice.

CHIN CONTROL SIMPLE			
ARA1609		A C S Chin Remote (Also select option ARA1601)	£113
ARA1289		Manual Swing Away Bracket RH	○
ARA1288		Manual Swing Away Bracket LH	○
DUAL CONTROL + BRACKET			
ARA1291		Dual Control RH (Please also select user control and bracket) (NC 1602, 1616, 1618, 1643, 1644, 1607, 1608, 1609)	£82
ARA1292		Dual Control LH (Please also select user control and bracket) (NC 1602, 1616, 1618, 1643, 1644, 1607, 1608, 1609)	£82
INVACARE IDC DUAL CONTROL + BRACKET			
ARA1645		Invacare IDC Dual control (NC1602, 1616)	£920
ARA1207		Standard Bracket RH NC (1609+1289)	○
ARA1208		Standard Bracket LH NC (1609+1288)	○
ARA1228		high and distance adjustable Bracket RH NC (1609+1289)	£41
ARA1229		high and distance adjustable Bracket LH NC (1609+1289)	£41
DRIVE UNIT OPTIONS			
FRAME COLOUR			
ARA0212	X	Flame Electric Red	◆
CENTRE WHEELS (DRIVING WHEELS)			
ARA0750	X	3.00 by 8.00 Grey Puncture Proof (Solid Inlay)	◆
PERFORMANCE OPTIONS			
MOTORTYPE			
ARA1284	X	8 km/h (4 pole)	◆
BATTERY			
ARA1260	X	73 Ah	◆
PERFORMANCE ACCESSORIES			
G-TRACK			
ARA1215		G-Track	£575
CHARGER			
ARA1343	X	Charger 8A (200 - 250V)	◆
SAFETY OPTIONS			
LIGHT			
ARA1367	X	Light system with indicators	◆
SAFETY ACCESSORIES			
TRANSPORTATION KIT			
ARA1494	X	4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19:2001) (Adequate headrest is compulsory, see options below)	◆
HEADRESTS ADEQUATE FOR PERSON TRANSPORTATION			
		INVACARE Headrest (Please select a headrest from this form)	○
ARA1495		Equivalent Headrest supporting the requirements of ISO7176-19:2001, provided independently from this order	○
FIXATION BELT			
ARA1329		Fixation Belt with Auto Lock	◆
ARA1332		H Belt (NC 1610)	£51
REARVIEW MIRROR			
ARA1415		RH (NC 1707)	£26
ARA1339		LH (NC 1708)	£26
ACCESSORIES			
CANEHOLDER			
ARA1312		RH	£16
ARA1363		LH	£16
TRAY			
ARA1404		Swing Away RH (NC 1207,1707)	£236
ARA1405		Swing Away LH (NC 1208,1708)	£236
TOTAL PRICE			

Order No:

Customer Ref:

LPF0099Q
TDX SP SSubject to Invacares standard terms conditions.
Invacare reserve the right to alter product specification without prior notice.

Non Confirming Items Guide

ARA1014	Powered Tilt (up to 28°) (max. user weight 160 kg)
ARA1021	Seat Lifter with electrical Tilt (max user weight 160 kg) 30cm rise 28° tilt
ARA0953	Pneumatic Recliner RH
ARA0954	Pneumatic recliner LH
ARA0947	Electric Recliner
ARA0556	Hanger , centre mounted, angle settable
ARA0557	Hanger, centre mounted, length compensating, electric. driven
ARA0533	Hanger Vari A
ARA0532	Hanger ADM 0° - 80°
ARA0531	Hanger ADE 0° - 80°
ARA1639	REM Type B (up to 4 electric actuators and light functionality)
ARA1207	Standard Bracket RH
ARA1208	Standard Bracket LH
ARA1707	Swing Away + height adjust. RH High End
ARA1708	Swing Away + height adjust. LH High End
ARA1644	REM Type B (up to 4 electric actuators and light functionality)
ARA1618	REM 550 with LCD screen
ARA1602	ACS REM24SD
ARA1616	G90A
ARA1608	A C S Centre Control (REM 24SD)
ARA1607	A C S Tray Control c/w Acrylic Tray & Actuator Functionality - DX REM 41
ARA1609	A C S Chin Remote
ARA1415	RH
ARA1339	LH
ARA1404	Swing Away RH
ARA1405	Swing Away LH

Special instructions (For consideration by Invacare's Engineering Department).

Special Instructions (For consideration by Approved Contractor).