

Invacare Top End All Terrain Wheelchair - UK Price List and Order Form

For accuracy of ordering, please fax order to Invacare Customer Service on 01656 776220 or e-mail to: ordersuk@invacare.com

www.topendwheelchair.co.uk



Provider name: _____ Account number: _____ PO number: _____
 Phone number: _____ Fax number: _____
 Address: _____

Customer name: _____ Phone number: _____
 Address: _____
 Ship to: Provider Customer (Order will be shipped to provider unless otherwise noted)

DATE:


QUOTE

ORDER

- Please indicate your choice with a cross in this square
- Indicate that this option has no upcharge


Manual Sport Wheelchair
 January 2012

Item	Description	Comments and/or non compatibilities information	Price
<input type="checkbox"/> CRFAT	Top End All Terrain wheelchair	Aluminium wheelchair built of 1-1/8" 6061 T6 Aluminium	£1,959
SEATING OPTIONS			
Seat width			
<input type="checkbox"/>	Seat width 14" (approx. 356 mm)		<input type="radio"/>
<input type="checkbox"/>	Seat width 15" (approx. 381 mm)		<input type="radio"/>
<input type="checkbox"/>	Seat width 16" (approx. 406 mm)		<input type="radio"/>
<input type="checkbox"/>	Seat width 17" (approx. 432 mm)		<input type="radio"/>
<input type="checkbox"/>	Seat width 18" (approx. 457 mm)		<input type="radio"/>
<input type="checkbox"/> CRF105	Seat width 19" (approx. 483 mm)		£67
<input type="checkbox"/> CRF105	Seat width 20" (approx. 508 mm)		£67
Seat depth			
<input type="checkbox"/>	Seat depth 14" (approx. 356 mm)		<input type="radio"/>
<input type="checkbox"/>	Seat depth 15" (approx. 381 mm)		<input type="radio"/>
<input type="checkbox"/>	Seat depth 16" (approx. 406 mm)		<input type="radio"/>
<input type="checkbox"/>	Seat depth 17" (approx. 432 mm)		<input type="radio"/>
<input type="checkbox"/>	Seat depth 18" (approx. 457 mm)		<input type="radio"/>
<input type="checkbox"/> CRF106	Seat depth 19" (approx. 483 mm)		£67
<input type="checkbox"/> CRF106	Seat depth 20" (approx. 508 mm)		£67
SEAT HEIGHTS OPTIONS			
Seat height front			
<input type="checkbox"/>	Seat height front 19" (approx. 483 mm)	Knee to heel range is 13.5"-17.5" for seat height front 21"	<input type="radio"/>
<input type="checkbox"/>	Seat height front 20" (approx. 508 mm)	Knee to heel range is 14.5"-18.5" for seat height front 21"	<input type="radio"/>
<input type="checkbox"/>	Seat height front 21" (approx. 533 mm)	Knee to heel range is 15.5"-19.5" for seat height front 21"	<input type="radio"/>
Seat height rear			
<input type="checkbox"/>	Seat height rear 15" (approx. 381 mm)	Only available with 24" wheels	<input type="radio"/>
<input type="checkbox"/>	Seat height rear 16" (approx. 406 mm)		<input type="radio"/>
<input type="checkbox"/>	Seat height rear 17" (approx. 432 mm)		<input type="radio"/>
<input type="checkbox"/>	Seat height rear 18" (approx. 457 mm)		<input type="radio"/>
<input type="checkbox"/>	Seat height rear 19" (approx. 483 mm)		<input type="radio"/>
<input type="checkbox"/>	Seat height rear 20" (approx. 508 mm)		<input type="radio"/>
CENTRE OF GRAVITY			
<input type="checkbox"/> CRF111	Active position 1.25" - 5.5"	Measured from the front end of the back post to the centre of the rear axle	<input type="radio"/>
<input type="checkbox"/> CRF113	Extended position -2.25" - 2.0"		£120
CAMBER SYSTEM			
<input type="checkbox"/> CRF630	0° camber tube with single camber		<input type="radio"/>
<input type="checkbox"/> CRF633	3° camber tube with single camber		<input type="radio"/>
<input type="checkbox"/> CRF636 (std)	6° camber tube with single camber		<input type="radio"/>
<input type="checkbox"/> CRF639	9° camber tube with single camber		<input type="radio"/>
<input type="checkbox"/> CRF640	0° extra camber tube		£100
<input type="checkbox"/> CRF643	3° extra camber tube		£100
<input type="checkbox"/> CRF646	6° extra camber tube		£100
<input type="checkbox"/> CRF649	9° extra camber tube		£100

FOOT SPACE WIDTH				
Standard	As per diagram		○	
				
FOOTPLATE OPTIONS				
<input type="checkbox"/>	Standard	Adjustable height and angle	○	
<input type="checkbox"/>	Standard	Calf and footstrap	○	
<input type="checkbox"/>	CRF191	Adjustable height/angle footrest with 2" risers	£64	
BACKREST				
Back type				
<input type="checkbox"/>	LKDWN (std)	Lock-down, folding ht/angle adj. Back	Manufactured of 1" aluminium tubing, only available with adjustable height (3/4" tube telescoping in 1" tube) ○	
<input type="checkbox"/>	CRF1070	Fixed back with fixed angle welded at 0°, adj ht	Manufactured of 1" aluminium available in fixed or adj height (7/8" telescoping in 1" tube), max height on fixed back is 16".	
<input type="checkbox"/>	CRF1073	Fixed back with fixed angle welded at 3°, adj ht		£133
<input type="checkbox"/>	CRF1076	Fixed back with fixed angle welded at 6°, adj ht		£133
<input type="checkbox"/>	CRF1079	Fixed back with fixed angle welded at 9°, adj ht		£133
Back height				
<input type="checkbox"/>	BH811	9"-11" (approx. 229 - 279 mm)	○	
<input type="checkbox"/>	BH1014	10"-14" (approx. 254 - 356 mm)	○	
<input type="checkbox"/>	BH1216	12"-16" (approx. 305 - 406 mm)	○	
<input type="checkbox"/>	BH1418	14"-18" (approx. 356 - 457 mm)	○	
<input type="checkbox"/>	BH1620	16"-20" (approx. 406 - 508 mm)	○	
Back upholstery and accessories				
<input type="checkbox"/>	CRFBAT	Ballistic nylon tension adj with straps and anchors	○	
<input type="checkbox"/>	CRF16	Fold over back upholstery	N/C with integral (CRF97) push handles and camouflage uphol. £67	
<input type="checkbox"/>	CRF30	Back pack	£32	
<input type="checkbox"/>	CRF97	Bolt-through integral push handles	N/C with fold over upholstery £68	
WHEELS				
Rear wheel options				
<input type="checkbox"/>	CRF24HPS	24" Sun® High Performance spoke	Supplied with 24" x 2" grey knobby tyres ○	
<input type="checkbox"/>	CRF25HPS	25" Sun® High Performance spoke	Supplied with 25" x 2.5" black knobby tyres ○	
<input type="checkbox"/>	CRF6724	24" Spinergy SPOX wheels	Supplied with 24" x 2" grey knobby tyres £328	
<input type="checkbox"/>	CRF6725	25" Spinergy SPOX wheels	Supplied with 25" x 2.5" black knobby tyres £328	
Spoke colour (only for Spinergy® wheels)				
<input type="checkbox"/>	SPCBK - Std	Black	If colour is not chosen, spokes will be black ○	
<input type="checkbox"/>	SPCY	Yellow	○	
<input type="checkbox"/>	SPCBW	White	○	
<input type="checkbox"/>	SPCR	Red	○	
<input type="checkbox"/>	SPCBL	Blue	○	
<input type="checkbox"/>	CRF27	Quad release axles	£84	
<input type="checkbox"/>	CRF76	Titanium axles	£102	
HANDRIMS				
Handrim options				
<input type="checkbox"/>	CRF77 (std)	Aluminium anodized handrims	○	
<input type="checkbox"/>	CRF75	Titanium handrims	£260	
<input type="checkbox"/>	CRF8	Coated handrims	£71	
<input type="checkbox"/>	CRF99	Natural Fit® standard grip handrims	£187	
<input type="checkbox"/>	CRF98	Natural Fit® super grip handrims	£217	
SEATING OPTIONS				
Seat accessories				
<input type="checkbox"/>	CRF29	Seat pouch	£23	
<input type="checkbox"/>	TED144	Heavy Duty seat pouch	£67	
<input type="checkbox"/>	CRF72	Positioning strap - Hook & Loop	£26	
<input type="checkbox"/>	CRF74	Positioning strap - Auto style	£33	
Cushions				
<input type="checkbox"/>	TED14	2" soft cushion	Filled with 2" soft foam which offer minimum protection £50	
Seat upholstery				
<input type="checkbox"/>	CRFATS	Adj tension seat with straps & anchor loop (standard)	○	

ARM OPTIONS			
Arm accessories			
<input type="checkbox"/>	CRF121	Fold-in style black aluminium side guards	Fold-in style flip back for easier transfers and are adjustable in height. Compatible with both high (push or pull to lock) or undermount (scissor) wheel locks. NC with adj height push handles if back height is 10-14" or less. £84
<input type="checkbox"/>	CRF19	Removable side guards without hole	£84
<input type="checkbox"/>	CRF18	Removable side guards with hole (quad)	£84
<input type="checkbox"/>	CRF20E	Fixed aluminium side guards even with wheels	Fixed back only £122
<input type="checkbox"/>	CRF201	Fixed aluminium side guards 1" above wheels	Fixed back only £122
<input type="checkbox"/>	CRF83E	Fixed carbon fibre side guards even with wheels	Fixed back only, N/C with tapered seat width £170
<input type="checkbox"/>	CRF831	Fixed carbon fibre side guards 1" above wheels	Fixed back only, N/C with tapered seat width £170
WHEEL LOCKS (choose one)			
<input type="checkbox"/>	TRA2	Push to lock wheel locks hi-mount	£71
MISCELLANEOUS ACCESSORIES			
<input type="checkbox"/>	CRF23	Rear anti tippers (E0971)	£64
<input type="checkbox"/>	CRF94	Alignment gauge	£47
<input type="checkbox"/>	RACE22	Top End T-Shirt, size L	£9
<input type="checkbox"/>	RACE22	Top End T-Shirt, size XL	£9
UPHOLSTERY COLOUR			
	Standard	Black	○
<input type="checkbox"/>	UPHI	Woodland camouflage	N/C with fold over back upholstery. Ballistic nylon tension adjustable with anchors and straps for the back upholstery. If a cushion is ordered it will also be camouflage. £67
<input type="checkbox"/>	UPH2	Desert camouflage	N/C with fold over back upholstery. Ballistic nylon tension adjustable with anchors and straps for the back upholstery. If a cushion is ordered it will also be camouflage. £88
FRAME COLOUR			
Contemporary colours			
<input type="checkbox"/>		Wet black	○
<input type="checkbox"/>		Black prism	○
<input type="checkbox"/>		White snow	○
<input type="checkbox"/>		Silver metallic	○
<input type="checkbox"/>		Stealth black	○
Way cool colours			
<input type="checkbox"/>		Electric blue	○
<input type="checkbox"/>		Midnight blue	○
<input type="checkbox"/>		Grasshopper green	○
Too hot colours			
<input type="checkbox"/>		Tangerine red	○
<input type="checkbox"/>		Electric red	○
<input type="checkbox"/>		Grape madness	○
<input type="checkbox"/>		Safety yellow	○
<input type="checkbox"/>		Sunkissed orange	○
<input type="checkbox"/>		Bubblegum pink	○

ADDITIONAL INFORMATION & MEASURING GUIDE



Top End Crossfire All Terrain Specifications

- Weight limit 250 pounds/113Kg
- Tapered Front End
- Frame Angle= 90 degrees
- Frame length= 25" for 14", 15" and 16" seat depths and seat depth plus 8" depths of 17"-20"
- Front wheels: 8" x 2" pneumatic with aluminium hubs
- Rear wheels: 25" x 2.5" black knobby pneumatic or 24" x 2" grey knobby pneumatic
- Castor centre to centre=footrest width plus 11"

Invacare Ltd: www.invacare.co.uk

Pencoed Technology Park, Pencoed, Bridgend, CF35 5AQ T: 01656 776222 F: 01656 776220

All rights reserved. Trademarks are identified by the symbols ™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted.

Sun® is a registered trademark of Sun Components Inc. Natural Fit is a registered trademark of Three Rivers Holding, LLC. Frog Legs is a registered trademark of Frog Legs, Inc.

Spinerig is a registered trademark for Spinerig Corp. Schwalbe is a registered trademark of RALF BOHLE GmbH in Reichshof, Germany.

Kenda is a registered trademark of Kenda Rubber Industrial Co. Ltd. Specifications and prices are subject to change without notification.

©2011 Invacare Corporation - Form No. 08-008

