



Yes, you can.®

e-spare cat Now Available Online
at www.invacare.co.uk

Account Name	Account No :
Address :	Order No :
Tel No:	Order Qty :
Delivery address :	Order Date :

Küschall® K - Junior
DEALER
Maximum user weight 60kg (9.5 Stone)

<p>Standard Features</p> <ul style="list-style-type: none"> Active rigid wheelchair Highly configurable Tension adjustable backrest upholstery Kuschall Starec Q/R Rear wheels Clothes guards 7 Different Seat Widths 7 Different Seat Depths Standard Base Price £1,927 	
--	--

Please make your selection with a cross in the box ◆ Standard feature ○ Free of Charge

Customer weight =		Kg	
Code	X	Description	Dealer
DDA0063		Küschall® K-series 2006	£1,927
DDA0000		Special construction # Please fax or email the full requirements #	P.O.A.

SEATING CHOICE

SEAT WIDTH (SB) *Air pump not mounted on 240-340mm widths*

DDA0114		SW 240 mm (NC: 0276)	○
DDA0116		SW 260 mm	○
DDA0118		SW 280 mm	○
DDA0120		SW 300 mm	○
DDA0122		SW 320 mm	○
DDA0124		SW 340 mm	○
DDA0126		SW 360 mm	○

SEAT DEPTH (ST)

DDA0401		Seat depth 250 mm	○
DDA0440		Seat depth 275 mm	○
DDA0406		Seat depth 300 mm	○
DDA0442		Seat depth 325 mm	◆
DDA0411		Seat depth 350 mm	○
DDA0414		Seat depth 375 mm	○
DDA0417		Seat depth 400 mm	○

BACK HEIGHT (RH)

In combination with armrests, back heights 270-330mm are only possible with standard push handles, buggy handle or without push handles

DDA0548		RH 270 mm (NC:1401)	○
DDA0580		RH 285 mm (NC:1401)	○
DDA0550		RH 300 mm (NC:1401)	○
DDA0503		RH 315 mm (NC:1401)	○
DDA0506		RH 330 mm (NC:1401)	○
DDA0509		RH 345 mm	○
DDA0512		RH 360 mm	○
DDA0514		RH 375 mm	◆
DDA0516		RH 390 mm	○
DDA0518		RH 405 mm	○
DDA0520		RH 420 mm	○
DDA0522		RH 435 mm	○
DDA0524		RH 450 mm	○
DDA0526		RH 465 mm	○
DDA0528		RH 480 mm	○


FRONT SEAT-TO-FLOOR (SHv)			
DDA0840		RSF 380 mm	○
DDA0841		RSF 390 mm	○
DDA0842		RSF 400 mm	○
DDA0843		RSF 410 mm	◆
DDA0844		RSF 420 mm	○
DDA0845		RSF 430 mm	○
DDA0846		RSF 440 mm	○
DDA0847		RSF 450 mm	○
REAR SEAT-TO-FLOOR (SHh)			
DDA0938		RSF 330 mm	○
DDA0939		RSF 340 mm	○
DDA0940		RSF 350 mm	○
DDA0941		RSF 360 mm	○
DDA0942		RSF 370 mm	◆
DDA0943		RSF 380 mm	○
DDA0944		RSF 390 mm	○
DDA0945		RSF 400 mm	○
DDA0946		RSF 410 mm	○
DDA0947		RSF 420 mm	○
DDA0948		RSF 430 mm	○
DDA0949		RSF 440 mm	○
DDA0950		RSF 450 mm	○
DDA0951		RSF 460 mm	○
DDA0952		RSF 470 mm	○
KNEE-TO-HEEL LENGTH (UL) 125mm (5") front wheels not possible with UL > 320mm. 100mm (4") not possible with UL >340mm			
Footrest high-mounted			
DDA1656		UL 150 mm	○
DDA1657		UL 160 mm	○
DDA1658		UL 170 mm	○
DDA1659		UL 180 mm	○
DDA1660		UL 190 mm	○
DDA1661		UL 200 mm	○
DDA1662		UL 210 mm	○
DDA1663		UL 220 mm	○
DDA1664		UL 230 mm	○
DDA1665		UL 240 mm	○
DDA1666		UL 250 mm	○
DDA1667		UL 260 mm	○
DDA1668		UL 270 mm	○
DDA1669		UL 280 mm	○
DDA1670		UL 290 mm	○
Standard footrest			
DDA1671		UL 300 mm	○
DDA1672		UL 310 mm	○
DDA1673		UL 320 mm	○
DDA1674		UL 330 mm	◆
DDA1675		UL 340 mm	○
DDA1676		UL 350 mm	○
DDA1677		UL 360 mm	○
DDA1678		UL 370 mm	○
DDA1679		UL 380 mm	○
DDA1680		UL 390 mm	○
DDA1681		UL 400 mm	○
DDA1682		UL 410 mm	○
BACK UPHOLSTERY			
DDA0701		Nylon upholstery with tension adjustable back	◆
PUSH HANDLES			
DDA0750		Without push handles	◆
DDA0702		Standard push handles	£93
DDA0703		Height adjustable push handles, rear set (NC :1411)	£141
DDA0704		Foldable push handles	£141
DDA0732		Children buggy push handle (NC: 0248,0251,0276,0280)	£177
SEAT UPHOLSTERY			
DDA0610		Standard, seat depth adjustable +25mm	◆

FOOTREST			
FOOTREST			
DDA1571		Footplate angle adjustable - Teddy bear (NC :1513)	○
DDA1572		Footplate angle adjustable - Küschall smiley (NC :1513)	◆
DDA1570		Footplate angle adjustable - Carbon footplate (NC :1513)	£141
DDA1575		Rear set footplate angle adjustable - Teddy bear (NC :1513)	£59
DDA1576		Rear set footplate angle adjustable - Küschall smiley (NC :1513)	£59
DDA1577		Rear set footplate angle adjustable - Carbon footplate (NC :1513)	£128
DDA1510		Foot tube rigid	○
DDA1513		Foot tube cover	£12
SIDE GUARDS			
Clothes guards/Mud guards			
DDA1401		Carbon clothes guards	◆
DDA1408		Carbon mud guards removable Adds 20mm to the total width if with 0° camber or Vario-Ax	£147
ARMREST (Armrest in combination with back heights 270mm - 330mm only available with standard push handles or without push handles)			
DDA1411		Tubular armrest, height adjustable, removable, swing-away, comfort padded	£133
CUSHION			
DDA0600		Seat cushion	£47
FRAME			
FRAME COLOUR, STANDARD			
DDA0300		Glossy Black	○
DDA0314		Black Matt	○
DDA0302		Royal Blue	○
DDA0303		Royal Red	○
DDA0350		Purple Fire	○
DDA0354		Night Shimmer	○
DDA0315		Silver Grey	○
DDA0352		Emerald Green	○
DDA0355		Midnight Blue	○
DDA0356		Ivory White	○
DDA0357		Umber Glaze	○
DDA0358		Smokey Grey	○
DDA0359		Flashy yellow	○
COLOUR-POWER-SET			
CPS not available for 508mm (20") rear wheels			
DDA0330		Techno-anthracite (NC: 0900)	£158
DDA0331		Lemon (NC : 0900)	£158
DDA0332		Apricot (NC : 0900)	£158
DDA0334		B-blue (NC : 0900)	£158
DDA0329		Comet Red (NC : 0900)	£158
FRONT WHEEL 125mm (5") front wheels not possible with UL > 320mm & 100mm (4") not possible with UL >340mm			
DDA0800		75mm (3") Supersport	○
DDA0801		100mm (4") Everyday	○
DDA0803		100mm (4") Starec techno-anthracite	£45
DDA0804		100mm (4") Starec lemon	£45
DDA0805		100mm (4") Starec apricot	£45
DDA0807		100mm (4") Starec b-blue	£45
DDA0808		100mm (4") Starec comet red	£45
DDA0809		100mm (4") Starec black	£45
DDA0824		100mm (4") Skater Wheel	£45
DDA0810		125mm (5") Light Wheel	◆
DDA0863		125mm (5") Comfort Wheel black rubber	○
DDA0813		125mm (5") Starec techno anthracite	£45
DDA0814		125mm (5") Starec lemon	£45
DDA0815		125mm (5") Starec apricot	£45
DDA0817		125mm (5") Starec b-blue	£45
DDA0818		125mm (5") Starec comet red	£45
DDA0819		125mm (5") Starec black	£45
QUICK RELEASE ADAPTER			
DDA0870		Quick release adapter	£60

REAR WHEEL			
REAR WHEEL FIXATION			
<i>On rear seat heights of 380mm or smaller, rear wheel extension is not possible with armrest</i>			
DDA0250		Standard rear wheel position (NC: 0263,0264,0265,0266,0267,0251,0248,0280)	◆
DDA0248		Light rear wheel extension (NC: 0266,0267,1705)	
DDA0251		Rear wheel extension (NC: 1705)	£88
CAMBER AXLE TYPE 1st AXLE			
DDA0276		Vario-Ax, adjustable Camber 1/3/7/10° (NC : 1302)	£184
DDA0277		Standard axle	◆
POSITION OF 1st AXLE			
<i>Axle position must only be selected when DDA0248 or DDA0251, Rear wheel extensions are chosen.</i>			
<i>2nd axle only applies when choosing DDA0251, the 5 position extension. If choosing a 2nd axle, 1st position must be selected on the 1st axle and 5th position on the 2nd axle</i>			
DDA0263		1st position from the front	○
DDA0264		2nd position from the front	○
DDA0265		3rd position from the front	○
DDA0266		4th position from the front	○
DDA0267		5th position from the front	○
CAMBER ANGLE 1st AXLE			
DDA0270		0° camber	○
DDA0271		3° camber	◆
DDA0272		6° camber	○
2nd CAMBER AXLE TYPE (When selecting a 2nd axle, 1st axle will be in 1st position & 2nd axle will be in 5th position)			
DDA0252		Standard	£180
DDA0280		Vario-Ax	£184
CAMBER ANGLE 2nd AXLE			
<i>Max. 3° difference in camber from 1st axle. If only 1 of the axles is a Vario-Axle the 2nd axle is only possible as 0° or 3° camber</i>			
DDA0270		0° camber	○
DDA0271		3° camber	○
DDA0272		6° camber	○
QUAD RELEASE AXLE			
DDA0909		Tetra Clip	£60
REAR WHEEL SIZE			
DDA0900		508mm (20") Starec rear wheel, with crossed spokes and light tyre. (NC: 0911,0904,0908,0913,0917,1104,1106,1202,1204,1207,1210,1211)	○
DDA0901		559mm (22") rear wheels (NC: 0904,0908,0913,0917,1004,1006,1010,1104,1106,1202,1207,1210,1211)	◆
DDA0902		609mm (24") rear wheels	○
REAR WHEEL TYPE			
DDA0904		609mm (24") Spider (NC :0910, 0911,1104) No spokeguards available	£158
DDA0905		508mm (20"), 559mm (22"), 609mm (24") Kuschall Starec wheel	◆
DDA0908		609mm (24") High Performance (NC: 0910,0911,1004,1006,1007,1104) No spokeguards available	£177
DDA0913		609mm (24") Mountain Bike wheel with 24" x 1.75" specific tyre (NC: 0910,0911,1104,1106,1202,1207,1210,1211)	£159
DDA0917		609mm (24") Spinergy wheel (NC: 0910,0911,1104) No spokeguards available	£409
<i>Only specify for the Kuschall Starec wheels</i>			
DDA0910		Crossed spoked	◆
DDA0911		Radial spoked	£93
REAR WHEEL TYRE			
DDA1001		508mm (20"), 559mm (22"), 609mm (24") light wheel tyre 7 bar (110 PSI)	◆
DDA1002		559mm (22"), 609mm (24") profil tyre 7 bar (110 PSI)	○
DDA1004		609mm (24") solid tyre, profil, grey	£47
DDA1006		609mm (24") solid tyre (KIK), black	£61
DDA1007		559mm (22"), 609mm (24") solid, light wheel	£47
DDA1010		609mm (24") Marathon Evolution" tyre (puncture protected with kevlar layer) (NC: 1020)	£71
ANTI-FLAT			
DDA1020		Anti-flat tyre tube protection	£14
HANDRIM			
DDA1102		508mm (20"), 559mm (22"), 609mm (24") Supergripp, distance 30mm	£114
DDA1104		609mm (24") Steel, distance 30mm	£44
DDA1106		609mm (24") Titanium, distance 30mm	£124
DDA1113		508mm (20"), 559mm (22"), 609mm (24") Light Alu, distance 30mm	○
DDA1114		508mm (20"), 559mm (22"), 609mm (24") Light Alu, distance 40mm	○
SPOKE GUARDS			
DDA1200		508mm (20"), 559mm (22"), 609mm (24") transparent	£82
DDA1202		609mm (24") Kuschall	£82
DDA1204		559mm (22"), 609mm (24") Boa-Boa	£115
DDA1207		609mm (24") "Original" Carbon	£122
DDA1210		609mm (24") Switzerland	£115
DDA1211		609mm (24") Pirates	£115

SAFETY			
BRAKES			
DDA1300		Standard brake lever bent	◆
DDA1302		Performance brake	£53
SAFETY OPTIONS			
BRAKE EXTENSION LEVER			
DDA1310		Extension lever for standard brake	£31
ANTI TIPPER			
DDA1702		Anti-tipper tube right	£31
DDA1703		Anti-tipper tube left	£31
DDA1704		Anti-tipper tube both sides	£61
DDA1705		Active Anti-tipper	£88
SEAT BELT			
DDA1730		Seat belt	£50
ACCESSORIES			
TRANSIT WHEELS			
DDA1700		Transit Wheels 75mm (3")	£147
PUMP			
DDA1715		Air Pump	£9
Küschall WATCH Mounted in castor housing			
DDA1712		Küschall watch, right side	£73
DDA1713		Küschall watch, left side	£73
TOOL KIT			
DDA1717		Tool Kit	◆
PRICE			
Total Price			£
Non Conforming Items Guide			
DDA0248		Light rear wheel extension	
DDA0251		Rear wheel extension	
DDA0263		1st axle position	
DDA0264		2nd axle position	
DDA0265		3rd axle position	
DDA0266		4th axle position	
DDA0267		5th axle position	
DDA0276		Vario-Ax, adjustable Camber 1/3/7/10°	
DDA0280		Vario-Ax	
DDA0900		20" rear wheels	
DDA0904		24" Spider Wheel	
DDA0908		24" High performance wheel	
DDA0910		Crossed spoked	
DDA0911		Radial spoked	
DDA0913		24" Mountain bike wheel	
DDA0917		24" Spinergy wheel	
DDA1004		24" solid tyre, profil, grey	
DDA1006		24" solid tyre (KIK), black	
DDA1007		22", 24" solid, light wheel	
DDA1010		Marathon Evolution tyres	
DDA1020		Anti-flat tyre tybe protection	
DDA1104		24" Steel, distance 3 cm	
DDA1106		24" Titanium, distance 3 cm	
DDA1202		24" küsch spokeguard	
DDA1204		22", 24" Boa-Boa spokeguard	
DDA1207		24" "Original" Carbon spokeguard	
DDA1210		24" Switzerland spokeguard	
DDA1211		24" Pirates spokeguard	
DDA1302		Performance brake	
DDA1401		Carbon clothes guards	
DDA1411		Tubular armrest, height adjustable, removable, swing-away, comfort padded	
DDA1513		Foot tube cover	
DDA1705		Active Anti-tipper	

Customer UL must be measured taking away the cushion thickness



UL = ULpat - SK

Possible combinations			
Rear wheel size	20"	22"	24"
Front wheel size	3"-4"	3"-5"	3"-5"
Front seat-to-floor height	38-40cm	40-42cm	43-45cm
Rear seat-to-floor height	33-42cm	36-44cm	38-47cm

Footplate advice :

Rear set for SD 25 - 30 cm
Standard for SD 32,5 - 40 cm

CPS-Set : Starec hubs ; 22" and 24" Starec Rims ; Forks ; Castor housing caps ; Front seat module caps ; Back angle joints ; Quick release axle caps for Starec hubs and Spinergy wheels.

Order No :

Customer Ref :

LPF0018T
Kuschall K-Junior
5 of 5

Subject to Invacare's standard terms conditions.
Invacare reserve the right to alter product specification without prior notice.