

INVACARE® TOP END®

Crossfire



Yes, you can.®

INVACARE® TOP END® Crossfire



The *Invacare Top End Crossfire* is the ultimate high-performance rigid custom manual wheelchair. The open frame construction lends the **Crossfire** its sleek, minimalist allure. If performance is what you are looking for, this is the chair for you.

With no underframe, the **Top End Crossfire** is easy to get in or out of a car. There are very few nuts and bolts to loosen or adjust, which translates to optimum performance and less maintenance and repairs. The **Crossfire Series** offers custom frame options such as a tapered seat width and the ergonomic seat option. The tapered seat width is great for people who want the chair to taper to their body and virtually disappear. The ergo seat option can improve posture, pressure distribution and propulsion for those who use 2-4 inches of seat slope. These wheelchairs are the perfect fit for anyone who wants a custom and light weight chair that is built for performance.

You may have noticed that disabled athletes never compete in folding wheelchairs. This is because rigid wheelchairs perform better at high levels of competition. Athletes and everyday people alike constantly seek to optimize their performance; a well-designed rigid wheelchair becomes part of the body, allowing easier access and freedom of movement. Body fit, aesthetics, independence, lightweight and durability are the main features of a rigid wheelchair, all of which results in a superior performance.

OPTIONS AND FEATURES



Each custom frame is built to accommodate a short front fork to reduce overall castor trail and keep the centre-of-gravity low to the ground for improved balance.



Two pressed bearings located in the fork allow for a more controlled rotation and help to eliminate flutter.



The side frame offers 4.25" (approx. 108 mm) of rear wheel centre-of-gravity adjustment.



Ergonomic seat option built right into the frame for better seat positioning.



TECHNICAL DATA



12" - 20"
(approx. 305 - 508 mm)



14" - 20"
(approx. 356 - 508 mm)



Front: 17" - 21"
(approx. 432 - 533 mm)
Rear: 15" - 21"
(approx. 381 - 533 mm)



Adj.: 9" - 20"
(approx. 229 - 508 mm)
Fixed: 8" - 16"
(approx. 203 - 406 mm)



Max. 250 lbs
(approx. 113 kg)

Invacare Limited
Pencoed Technology Park - Pencoed
Bridgend CF35 5AQ - United Kingdom
Tel: +44 1 656 776 222
Fax: +44 1656 776 220
E-mail: uk@invacare.com
Sales Order E-mail:
ordersuk@invacare.com
www.invacare.co.uk

Invacare Ireland Ltd
Unit 5 - Seatown Business Campus
Seatown Road - SWORDS
County Dublin - Ireland
Tel. +353 1 8107084
Fax +353 1 8107085
Email: Ireland@invacare.com
www.invacare.ie

© 2012 Invacare International Sarl
All rights reserved.

All information quoted is believed
to be correct at time of print.
Invacare reserves the right to alter product
specifications without prior consultation.

TopEnd Crossfire - UK - 03/2012 - 1524386

Please refer to order form for all
technical specifications.



Yes, you can.®