



küschall[®] *Compact*

Active wheelchair

New!

Precision is essential for *küschall*, and pushes us to bring more and more information to you. This is the reason why this new portfolio edition includes the weight of each part.

Thanks to these details, you will be able to compare the weight of all options and make a safe choice; you will then easily estimate the whole weight of the *Compact* wheelchair you have selected.

Except any other indication, all weights indicated are based on seat width 40cm and seat depth 40cm. Options weights are calculated by pair when requested (e.g. for the wheels, push handles, brakes...) and with their related fixations.

Compact

Contents

- Frame
- Backrest.....
- Seating
- Footrest
- Side guards
- Cushions and options
- Colours
- Front wheels
- Rear wheels.....
- Brakes
- Safety options.....
- Accessories

Frame material and angle

DDP0211

Standard frame






DDP0224
DDP0205
DDP0219

Abducted frame,
left, right or both sides






Backrest

Backrest upholstery

DDP0700	Standard, light backrest upholstery	
	<i>320 grams</i>	
DDP0701	Black backrest, tension adjustable with <i>Velcro®</i>	
	<i>560 grams (including Velcro® straps)</i>	
DDP0735	Aero backrest, tension adjustable with <i>Velcro®</i>	
	<i>600 grams (including Velcro® straps)</i>	

Back tube

DDP0705	Lumbal tubes	
DDP0734	Foldable back (to the front), angle adjustable	
	<i>1'080 grams (including push handles)</i>	
DDP0707 DDP0718 DDP0709 DDP0719 DDP0711 DDP0712	Back angle adjustable Adjusted: 82°, 86°, 90°, 94°, 98°, 102°	

Push handles

DDP0750	Without push handle	
DDP0702	Push handles fixed, long	
	<i>240 grams the pair (without back tube)</i>	
DDP0751	Push handles fixed, mini	
	<i>160 grams the pair (without back tube)</i>	
DDP0703	Push handles, height adjustable, rear set	
	<i>1'800 grams the pair (incl. tubes and fixations)</i>	
DDP0733	Push handles, height adjustable, integrated	NO PICTURE AVAILABLE
	<i>980 grams the pair (incl. tubes and fixations)</i>	
DDP0704	Push handles, foldable	
	<i>1'120 grams the pair (incl. tube and fixations)</i>	

Backrest option

DDP1491

Backrest stabilization bar,
swing away and removable

220 grams




Seating

Seat upholstery


DDP0611 Tension adjustable

560 grams (including Velcro® bands and fixations)




DDP0613 Seat upholstery with central Velcro® band
(soft band)

560 grams (including Velcro® bands and fixations)



DDP0614 Seat upholstery with lateral Velcro® bands
(rough bands)

560 grams (including Velcro® bands and fixations)



Footrest

Footrests

DDP1562
DDP1563


Two pieces footrest, plastic,
set to the front or to the back



580 grams the pair (incl. fixations)

DDP1564
DDP1565

Two pieces footrest, angle adjustable,
set to the front or to the back



900 grams the pair (incl. fixations)

DDP1578
DDP1579


Two pieces carbon footplate,
set to the front or to the back



740 grams the pair (incl. fixations)

DDP1566
DDP1567

One piece footrest flip-up,
angle adjustable, fixed on one side when
folded,
set to the back or to the front



850 grams (incl. fixations)

Legrests

DDP1540 Legrest,
angle open, 70°



720 grams the pair

DDP1541 Legrest,
angle open, 80°



720 grams the pair

DDP1542 Legrest,
angle open, 90°



720 grams the pair

Foot and leg supports

DDP1534 Calf-strap



30 grams

DDP1530 Heel loops
(only with 2 pieces plastic footplates)



90 grams the pair




DDP1536 Instep loops
(only with 2 piece plastic footplates)






60 grams the pair

Side guards

Clothes guards and mud guards

DDP1400	Carbon clothes-guards, fixed	
	<i>Short: 520 grams the pair (incl. fixations)</i> <i>Large: 680 grams the pair (incl. fixations)</i>	
DDP1408	Carbon mud guards, removable (This option adds 2 cm to the total width of the wheelchair)	
	<i>Short: 660 grams the pair (incl. fixations)</i> <i>Large: 824 grams the pair (incl. fixations)</i>	
DDP1405	Carbon mud guards, fixed	
	<i>Short: 460 grams the pair (incl. fixations)</i> <i>Large: 620 grams the pair (incl. fixations)</i>	

Siderests

DDP1416 DDP1421	Siderest height adjustable (with tools), removable	
	<i>Short: 1'580 grams the pair</i> <i>Long: 1'700 grams the pair</i>	
DDP1430 DDP1431	Siderest flip-back height adjustable, removable, in frame color (This option adds 2 cm to the total width of the wheelchair)	
	<i>Short: 2'340 grams the pair</i> <i>Long: 2'480 grams the pair</i>	
DDP1440 DDP1441	Siderest flip-back, height adjustable, removable, with desk function, (This option adds 2 cm to the total width of the wheelchair)	
	<i>Short: 2'180 grams the pair</i> <i>Long: 2'760 grams the pair</i>	

Armrests

DDU1411

Tubular armrest, height adjustable, removable, swing away, comfort padded



820 grams the pair

Cushions and options

DDP0660 Flo-tech Lite cushion



*1'020 grams (SB40 cm*ST42 cm)*

DDP0600 Seat cushion

NO PICTURE AVAILABLE

*510 grams (SB40 cm*ST40 cm)*

DDP0606 Seat bag, click in and removable








130 grams

Colours

Frame colours

Glossy black		Black matt	
Royal blue		Royal red	
Purple fire		Silver grey	
Emerald green		Midnight blue	
Umber glaze		Ivory white	
Smokey grey		Night shimmer	

CPS-colours

Techno-anthracite		Lemon	
B-blue		Comet red	
Apricot		Copper glaze	
Artic blue			

The Colour Power Set (CPS) applies on the following parts and options: Starec & High Performance rear wheel hubs with quick release protection, 22" and 24" Starec rear wheel rims, Starec front wheels, forks, front seat module caps, rear wheel adapter plate.

Front wheels

DDP0800	Supersport	
	<i>3": 212 grams the pair</i>	
	Starec wheel, in various CPS colours	
	<i>4": 364 grams the pair</i> <i>5": 536 grams the pair</i>	
DDP0863	Comfort wheel with black rubber	
	<i>5": 780 grams the pair</i>	
DDP0810	Light wheel	
	<i>5": 624 grams the pair</i> <i>6": 788 grams the pair</i>	
DDP0824	Skater wheel	
	<i>4": 344 grams the pair</i>	
DDP0801	Everyday wheel	
	<i>4": 548 grams the pair</i>	
DDP0821	Pneumatic wheel	
	<i>6": 616 grams the pair</i>	
DDP0811	Softroll	
	<i>5": 860grams the pair</i>	

Rear wheels

Rear wheel fixation and position

DDP1118
DDP1119

One arm drive system,
pneumatic tyres and alu handrims,
right or left side


5'630 grams (incl. wheels, tyres and handrims)



DDP0250

Rear wheel standard fixation
Rear wheel positions


240 grams the pair



DDP0249

Rear wheel extension
Rear wheel positions

760 grams the pair



DDP0909

Quad or Tetra quick release axle








440 grams the pair



Rear wheel type

DDP0915	küschall Starec wheel <i>22": 1'456 grams the pair</i> <i>24": 1'592 grams the pair</i> <i>26": 1'700 grams the pair</i>	
DDP0908 DDP0918	High Performance wheel (standard hub is black – if no CPS choice) <i>24": 1'356 grams the pair</i> <i>25": 1'372 grams the pair</i>	
DDP0913	Mountain Bike wheel <i>24": 1'712 grams the pair</i>	
DDP0917	Spinergy wheel <i>24": 1'476 grams the pair</i> <i>25": 1'508 grams the pair</i>	
DDP0904	Spider wheel <i>24": 2'168 grams the pair</i>	

Rear wheels tyres

<p>DDP1001 Light wheel 7 bar</p> <p>22": 736 grams the pair 24": 816 grams the pair 26": 844 grams the pair</p>		<p>DDP1002 Profile tyre 7 bar</p> <p>22": 592 grams the pair 24": 624 grams the pair</p>	
<p>DDP1007 Solid, light wheel</p> <p>22": 1'032 grams the pair 24": 1'032 grams the pair</p>	<p>NO PICTURE AVAILABLE</p>	<p>DDP1004 Solid tyre, profile, grey</p> <p>24": 1'140 grams the pair</p>	
<p>DDP1003 Sports hall</p> <p>24": 640 grams the pair 26": 728 grams the pair</p>		<p>DDP1006 Solid tyre (KIK Type), black</p>	
<p>DDP1010 Marathon Evolution Plus</p> <p>24": 1'024 grams the pair 25": 1'076 grams the pair 26": 1'100 grams the pair</p>		<p>DDP1011 Schwalbe speedrun</p> <p>24": 356 grams the pair 25": 996 grams the pair</p>	

Handrims

<p>DDP1102 Supergripp</p> <p>22": 1'060 grams the pair 24": 1'336 grams the pair 26": 1'376 grams the pair</p>		<p>DDP1104 Steel</p> <p>24": 1'276 grams the pair 26": 1'392 grams the pair</p>	
<p>DDP1106 Titanium</p> <p>24": 808 grams the pair</p>		<p>DDP1109 Maxgripp</p> <p>24": 1'160 grams the pair</p>	
<p>DDP1110 Tetra</p> <p>24": 2'092 grams the pair 26": 2'284 grams the pair</p>		<p>DDP1113, DDP1114 Light Alu</p> <p>22": 592 grams the pair 24": 716 grams the pair 25": 744 grams the pair 26": 760 grams the pair</p>	

Spoke guards

DDP1200	Transparent	
	<i>22": 480 grams the pair 24": 660 grams the pair 26": 800 grams the pair</i>	
DDP1207	«Original» Carbon	
	<i>24": 720 grams the pair</i>	
DDP1211	Pirates	
	<i>24": 660 grams the pair</i>	
DDP1202	küsch	
	<i>24": 680 grams the pair</i>	
DDP1204	Boa-Boa	
	<i>22": 520 grams the pair 24": 680 grams the pair</i>	
DDP1210	Switzerland	
	<i>24": 660 grams the pair</i>	

Brakes

DDP1300	Standard push brakes, levers bent	
	<i>380 grams the pair</i>	
DDP1301	Push brakes, levers straight	
	<i>380 grams the pair</i>	
DDP1302	Performance brakes	
	<i>420 grams the pair</i>	
DDP1303	Active brakes	
	<i>300 grams the pair</i>	
DDP1306	Pull brakes	
	<i>400 grams the pair</i>	
DDP1304	Attendant drum brakes	NO PICTURE AVAILABLE
DDP1310	Extension lever for standard brakes	
	<i>Same weight as without extension lever</i>	

Safety options

DDF1730 Seat belt



230 grams

DDF1710 Passive lights
DDF1711

NO PICTURE AVAILABLE

60 grams the pair

DDF1702 Anti tipper, height adjustable,
DDF1703 active and inactive position, removable,
DDF1704 right, left or both sides








Short: 370 grams (without fixations)
Long: 400 grams (without fixations)

Fixations:

210 grams with rear wheel extension

Accessories

DDP1706 DDP1707	Cane-holder right or left	
	<i>240 grams</i>	
DDP1708 DDP1709	Tipper-Aid right or left	
	<i>180 grams</i>	
DDP1700	Transit wheels 3"	
	<i>1'160 grams (including fixations)</i>	
DDP1717	Tool kit	
	<i>60 grams</i>	
DDP1715	Air pump	
	<i>100 grams</i>	

