



Yes, you can:

e-spare cat Now Available Online at
www.invacare.co.uk

Customer Name :	Account No :
Dealer :	Order No :
Address :	Order Qty :
Tel No:	Order Date :
Delivery address :	

Küschall® Ultra-Light
DEALER
Maximum user weight 130kg (20.5 Stone)

Standard Features

Lightweight aluminium folding frame
Highly configurable positioning
Clothes guards
Starec Rear Wheels
Aluminium handrims
2 Piece footplates
12 Different seat width options
10 Different seat depth options
Standard Base Price £1,307

Please make your selection with a cross in the box ◆ Standard feature ○ Free of Charge

Customer weight =		Kg	Dealer
Code	X	Description	
DDU0010		Küschall® Ultra-Light	£1,307
DDU0000		Special Construction # Please fax or email the full requirements #	P.O.A.

SEATING CHOICE

FRONT FRAME		Minimum seat angle (SHv-SHh) with dynamic frame is 30mm	
DDU0203		Active, 70°	◆
DDU0212		Dynamic, 80° (NC: 0820, 0821, 0822)	£12

Note: On adducted frames: Foot width = seat width less 140mm. Maximum is 290mm

Folding back function is limited when selected with Adducted frames		
DDU0201	Adducted Frame 70° (NC: 0118,0120)	○
DDU0243	Adducted Frame 80° (NC: 0118,0120)	○

SEAT WIDTH (SB) Note: Allow an extra 10mm on the actual seat width required.

Note: Seat widths 280mm - 320mm are only compatible with 1 or 2 piece carbon footplates

Note: On Adduction frames, Seat widths less than 380mm are only compatible with 1 or 2 piece carbon footplates

DDU0118	Seat width 280 mm	single cross brace Max user weight is - 130kg (20.5 Stone)	(NC: 0660)	○
DDU0120	Seat width 300 mm	single cross brace Max user weight is - 130kg (20.5 Stone)	(NC: 0660)	○
DDU0122	Seat width 320 mm	single cross brace Max user weight is - 130kg (20.5 Stone)	(NC: 0660)	○
DDU0124	Seat width 340 mm	single cross brace Max user weight is - 130kg (20.5 Stone)	(NC: 0660)	○
DDU0126	Seat width 360 mm	single cross brace Max user weight is - 130kg (20.5 Stone)	(NC: 0660)	○
DDU0128	Seat width 380 mm	single cross brace Max user weight is - 130kg (20.5 Stone)		○
DDU0130	Seat width 400 mm	single cross brace Max user weight is - 130kg (20.5 Stone)		○
DDU0132	Seat width 420 mm	single cross brace Max user weight is - 130kg (20.5 Stone)		○
DDU0134	Seat width 440 mm	double cross brace Max user weight is 130kg (20.5 Stone)	(NC: 0660)	○
DDU0136	Seat width 460 mm	double cross brace Max user weight is 130kg (20.5 Stone)	(NC: 0660)	○
DDU0138	Seat width 480 mm	double cross brace Max user weight is 130kg (20.5 Stone)	(NC: 0660)	○
DDU0140	Seat width 500 mm	double cross brace Max user weight is 130kg (20.5 Stone)	(NC: 0660)	○

SEAT DEPTH (ST)

DDU0408	Seat depth 320 mm	(NC: 0660)	○
DDU0410	Seat depth 340 mm	(NC: 0660)	○
DDU0412	Seat depth 360 mm	(NC: 0660)	○
DDU0415	Seat depth 380 mm		○
DDU0417	Seat depth 400 mm		◆
DDU0420	Seat depth 420 mm		○
DDU0422	Seat depth 440 mm		○
DDU0424	Seat depth 460 mm	(NC: 0660)	○
DDU0427	Seat depth 480 mm	(NC: 0660)	○
DDU0429	Seat depth 500 mm	(NC: 0660)	○

BACK HEIGHT (RH)				
DDU0550		RH 300 mm	Straight backposts supplied as standard (NC: 1491,0734)	○
DDU0503		RH 315 mm	Straight backposts supplied as standard (NC: 1491,0734)	○
DDU0506		RH 330 mm	Straight backposts supplied as standard (NC: 1491,0734)	○
DDU0509		RH 345 mm	(NC: 1491)	○
DDU0512		RH 360 mm	(NC: 1491)	○
DDU0514		RH 375 mm	(NC: 1491)	◆
DDU0516		RH 390 mm	(NC: 1491)	○
DDU0518		RH 405 mm		○
DDU0520		RH 420 mm		○
DDU0522		RH 435 mm		○
DDU0524		RH 450 mm		○
DDU0526		RH 465 mm		○
DDU0528		RH 480 mm		○
DDU0530		RH 495 mm		○
DDU0532		RH 510 mm		○
FRONT SEAT-TO-FLOOR HEIGHT (SHv) Max Front height with Dynamic frame is 520mm!				
DDU0845		FSF 430 mm		○
DDU0846		FSF 440 mm		○
DDU0847		FSF 450 mm		○
DDU0848		FSF 460 mm		○
DDU0849		FSF 470 mm		○
DDU0850		FSF 480 mm		○
DDU0851		FSF 490 mm		○
DDU0852		FSF 500 mm		◆
DDU0853		FSF 510 mm		○
DDU0854		FSF 520 mm		○
DDU0855		FSF 530 mm		○
REAR SEAT-TO-FLOOR HEIGHT (SHh) Please note: SHh can be maximum 100mm lower, and a maximum 10mm higher than SHv				
Mudguards in combination with a rear wheel extension are not compatible with rear seat to floor heights of 410mm or lower				
DDU0944		RSF 390 mm	(NC: 1405 + 0249)	○
DDU0945		RSF 400 mm	(NC: 1405 + 0249)	○
DDU0946		RSF 410 mm	(NC: 1405 + 0249)	○
DDU0947		RSF 420 mm		○
DDU0948		RSF 430 mm		○
DDU0949		RSF 440 mm		◆
DDU0950		RSF 450 mm		○
DDU0951		RSF 460 mm		○
DDU0952		RSF 470 mm		○
DDU0953		RSF 480 mm		○
DDU0954		RSF 490 mm		○
DDU0955		RSF 500 mm		○
KNEE-TO-HEEL LENGTH Please note that for knee to heel lengths less than 390mm and depending on the configuration the footplate may be upmounted. On 70° frame, knee to heel = (SHv - 20mm), on 80° frame, knee to heel = (SHv - 30mm)				
DDU1661		UL 200 mm	(NC: 1564,1565)	○
DDU1662		UL 210 mm	(NC: 1564,1565)	○
DDU1663		UL 220 mm	(NC: 1564,1565)	○
DDU1664		UL 230 mm	(NC: 1564,1565)	○
DDU1665		UL 240 mm	(NC: 1564,1565)	○
DDU1666		UL 250 mm	(NC: 1564,1565)	○
DDU1667		UL 260 mm	(NC: 1564,1565)	○
DDU1668		UL 270 mm	(NC: 1564,1565)	○
DDU1669		UL 280 mm		○
DDU1670		UL 290 mm		○
DDU1671		UL 300 mm		○
DDU1672		UL 310 mm		○
DDU1673		UL 320 mm		○
DDU1674		UL 330 mm		○
DDU1675		UL 340 mm		○
DDU1676		UL 350 mm		○
DDU1677		UL 360 mm		○
DDU1678		UL 370 mm		○
DDU1679		UL 380 mm		○
DDU1680		UL 390 mm		○
DDU1681		UL 400 mm		○
DDU1682		UL 410 mm		○

DDU1683		UL 420 mm	○
DDU1684		UL 430 mm	◆
DDU1685		UL 440 mm	○
DDU1686		UL 450 mm	○
DDU1687		UL 460 mm	○
DDU1688		UL 470 mm	○
DDU1689		UL 480 mm	○
DDU1690		UL 490 mm	○
DDU1691		UL 500 mm	○
DDU1692		UL 510 mm	○
BACKREST CHOICE			
DDU0700		Standard, light upholstery (NC: 0750, 1491, 0703)	◆
DDU0722		Tension velcro adjustable back, black & grey	£97
DDU0735		Aero back, tension velcro adjustable	£158
PUSH HANDLES			
DDU0702		Push-handles fixed, long	◆
DDU0751		Push handles fixed, mini (NC 1304)	○
DDU0750		Without push handles	○
DDU0703		Height adjustable, rear set	£180
DDU0733		Height adjustable, integrated	£87
DDU0704		Foldable push handles. (NC 1304,1491)	£140
BACK TUBE OPTION If no option selected, straight backrest tubes will be assembled			
DDU0705		Lumbal tubes (NC: 0734)	○
DDU0720		Angle adjustable back (NC: 0734) (Backrest angle must be selected from the 6 options below)	£143
DDU0734		Fold down back (folding forward onto the seat) (NC: 0712,1411) (Backrest angle must be selected, see table on back page for angle possibilities)	£52
BACKREST ANGLE			
DDU0707		Angle adjustable backrest, adjusted 82°	○
DDU0718		Angle adjustable backrest, adjusted 86°	○
DDU0709		Angle adjustable backrest, adjusted 90°	◆
DDU0719		Angle adjustable backrest, adjusted 94°	○
DDU0711		Angle adjustable backrest, adjusted 98°	○
DDU0712		Angle adjustable backrest, adjusted 102°	○
BACKBRACE			
DDU1491		Backbrace swing away and removable	£88
SEATING			
DDU0611		Standard adjustable seat upholstery with 1 central velcro band	◆
FOOTREST			
Seat width 280-320mm & Adduction frame with seat width less than 380mm only 1pc or 2pc Carbon footplates are available			
DDU1562		2 pc. plastic, set to front	◆
DDU1563		2 pc. plastic, set to back	○
DDU1564		2 pc. angle adjustable set to front	£92
DDU1565		2 pc. angle adjustable set back	£92
DDU1578		2 pc. carbon footplates, set to the front (NC: 1530,1536)	£237
DDU1579		2 pc. carbon footplates, set to the rear (NC: 1530,1536)	£237
DDU1566		1 pc. carbon Flip-up, adjustable, fixed on 1 side when folded, set to front (NC: 1530,1536)	£118
DDU1567		1 pc. carbon Flip-up, adjustable, fixed on 1 side when folded, set to back (NC: 1530,1536)	£118
DDU1539		Foot strap	£19
FOOT AND LEG SUPPORT			
DDU1533		Calf strap	◆
DDU1530		Heel loops	○
DDU1536		Instep loops	○
SIDE GUARDS			
CLOTHES & MUD GUARDS			
DDU1400		Clothes guards carbon, fixed (NC:1416,1421,1430,1431,1440,1441)	◆
DDU1405		Mud guards carbon, fixed (NC: 1416,1421,1430,1431,1440,1441)	£112
DDU1408		Mudguards carbon, removable (NC:1416,1421,1430,1431,1440,1441) (Adds 20mm to the total width of the chair)	£199
ARMRESTS (Includes clothes guards) (Adds 20mm to total width of the chair)			
DDU1416		Armrest height adjustable (with tools), removable, short	○
DDU1421		Armrest height adjustable (with tools), removable, long	○
DDU1430		Armrest flip-back height adjustable removable, short (in frame colour, with carbon clothes guards)	£293
DDU1431		Armrest flip-back height adjustable removable, long (in frame colour, with carbon clothes guards)	£293
DDU1440		Armrest with desk function, flip-back, height adjustable, removable, short (NC: 1300,0913,0283)	£150
DDU1441		Armrest with desk function, flip-back, height adjustable, removable, long (NC: 1300,0913,0283)	£180
DDU1411		Tubular armrests height adjustable,removable,swingaway, comfort padded armpad, no clothesguards	£170
If Tubular armrests are selected, only 1 other adaption can be fitted. Adaptions include Anti tips, Tipper tube, Cane holder and Transit wheels			

SEATING OPTIONS			
SEAT CUSHION			
DDU0600		Standard seat cushion	£47
DDU0660		Flo Tech Lite seat cushion. Only available for seat width 380 to 420mm and for seat depth 380 to 440mm	£101
DDU0606		Seat Bag (Click in and removable)	£36
FRAME COLOUR			
DDU0300		Glossy-Black	○
DDU0314		Black Matt	○
DDU0356		Ivory White	○
DDU0302		Royal Blue	○
DDU0303		Royal Red	○
DDU0350		Purple Fire	○
DDU0354		Night Shimmer	○
DDU0315		Silver Grey	○
DDU0352		Emerald Green	○
DDU0355		Midnight Blue	○
DDU0357		Umber Glaze	○
DDU0358		Smokey Grey	○
DDU0359		Flashy Yellow	○
DDU0363		Black-Grey Rustico	○
COLOUR-POWER-SET			
DDU0330		Colour-Power-Set techno-anthracite	£158
DDU0331		Colour-Power-Set lemon	£158
DDU0332		Colour-Power-Set apricot	£158
DDU0334		Colour-Power-Set b-blue	£158
DDU0329		Colour-Power-Set comet red	£158
DDU0337		Colour-Power-Set copper glaze	£158
DDU0336		Colour-Power-Set artic blue	£158
DDU0338		Colour-Power-Set chrome	£158
FRONT WHEEL ATTACHMENT			
DDU0864		High mounted (Frame will be 30mm lower)	○
DDU0865		Standard mounted	◆
DDU0866		Low mounted (Frame will be 30mm higher)	○
FRONT WHEEL 6" & 7" not possible with Dynamic 80° frame			
DDU0800		75mm (3") Supersport	○
DDU0801		100mm (4") Everyday	£45
DDU0803		100mm (4") Starec techno-anthracite	£45
DDU0804		100mm (4") Starec lemon	£45
DDU0805		100mm (4") Starec apricot	£45
DDU0807		100mm (4") Starec b-blue	£45
DDU0808		100mm (4") Starec comet red	£45
DDU0830		100mm (4") Starec copper glaze	£45
DDU0831		100mm (4") Starec artic blue	£45
DDU0809		100mm (4") Starec black	£45
DDU0861		100mm (4") Starec chrome	£45
DDU0824		100mm (4") Skater Wheel	£45
DDU0810		125mm (5") Light Wheel	○
DDU0863		125mm (5") Comfort Wheel black rubber	◆
DDU0813		125mm (5") Starec techno anthracite	£45
DDU0814		125mm (5") Starec lemon	£45
DDU0815		125mm (5") Starec apricot	£45
DDU0817		125mm (5") Starec b-blue	£45
DDU0818		125mm (5") Starec comet red	£45
DDU0833		125mm (5") Starec copper glaze	£45
DDU0832		125mm (5") Starec artic blue	£45
DDU0819		125mm (5") Starec black	£45
DDU0862		125mm (5") Starec chrome	£45
DDU0820		150mm (6") Light wheel	○
DDU0821		150mm (6") Pneumatic Wheel	£18
DDU0822		175mm (7") Light Wheel	○
REAR WHEEL ADAPTION PLATE			
DDU0250		Standard fixation (NC: 0267, 0268)	◆
DDU0249		Rear wheel extension (Black in colour)	£47

REAR WHEEL POSITION			
DDU0263		1st position (NC: 1304)	○
DDU0264		2nd position	○
DDU0265		3rd position	◆
DDU0266		4th position	○
DDU0267		5th position	○
DDU0268		6th position	○
QUAD QUICK RELEASE AXLE			
DDU0909		Quick release axle (NC: 1304)	£60
REAR WHEEL SIZE			
DDU0901		559mm (22") (NC: 0904,0908,0918,0913,0917,0916)	○
DDU0902		609mm (24")	◆
DDU0925		635mm (25") (NC: 0916,0904,0915,0908,0913,0283,1102,1104)	○
DDU0903		660mm (26") (NC: 0283,0916,0904,0908,0918,0913,0917)	○
REAR WHEEL CAMBER (3° Not available with drum brakes and performance brakes combined)			
DDU0282		Camber 0°	◆
DDU0283		Camber 3° (NC: 1302,1304,0916)	○
REAR WHEEL TYPE			
DDU0904		609mm (24") Spider-Wheel (NC: 0910,0911,1104,1120) No spokeguards available	£158
DDU0915		559mm (22"), 609mm (24"), 660mm (26") Kuschall Starec	◆
DDU0908		609mm (24") High Performance (Starec hub) (NC: 0910,0911,1004,1006,1007,1104) No spokeguards available	£177
DDU0918		635mm (25") High Performance complete wheel, fitted with Marathon tyre & Aluminium handrim (Starec hub) No spokeguards available (NC: 0910,0911)	£258
DDU0913		609mm (24") Mountain Bike with 24" x 1.75" specific tyre. (NC,0910,0911,1104,1106,1120) Available with Transparent or Boa-Boa spokeguards only	£159
DDU0917		609mm (24") Spinergy wheel (NC: 0910,0911) No spokeguards available on 24" or 25" Spinergys 635mm (25") Spinergy wheel (NC: 0910,0911,1001,1002,1004,1006,1007,1102,1104,1106)	£409
TYPES OF SPOKES (Select only with Kuschall Starec wheels)			
DDU0910		Cross spoked	○
DDU0911		Radial spoked	○
E-MOTION ADAPTATION (Does not include e-motion wheels) Only compatible with 609mm (24") wheels & 0 Degree camber			
DDU0916		E-motion adaptation plate (NC: 1702,1703,1706,1707)	£182
REAR WHEEL TYRE			
DDU1001		559mm (22"), 609mm (24"), 660mm (26") Light wheel 7 bar	◆
DDU1002		559mm (22"), 609mm (24") Profile tyre 7 bar	○
DDU1004		609mm (24") Solid tyre, profil, grey	£47
DDU1006		609mm (24") Solid tyre (KIK type), black	£61
DDU1007		559mm (22"), 609mm (24") Solid, light wheel	£47
DDU1010		609mm (24"), 635mm (25"), 660mm (26") Marathon Evolution tyre from Schwalbe (NC: 1020)	£71
DDU1011		609mm (24"), 635mm (25") Schwalbe speedrun	£60
ANTI-FLAT			
DDU1020		Anti Flat tyre tube protection	£14
HANDRIM			
DDU1113		559mm (22"), 609mm (24"), 635mm (25"), 660mm (26") Light Alu, dist. 3 cm	◆
DDU1114		559mm (22"), 609mm (24"), 635mm (25"), 660mm (26") Light Alu, dist. 4 cm	○
DDU1102		559mm (22"), 609mm (24"), 660mm (26") Supergrripp, dist. 3 cm	£114
DDU1104		609mm (24"), 660mm (26") Steel, dist. 3 cm	£44
DDU1106		609mm (24") Titanium, dist. 3 cm	£124
DDU1120		609mm (24"), 635mm (25"), 30mm Surge LT, 3cm distance	£230
SPOKE GUARDS			
DDU1200		559mm (22"), 609mm (24"), 660mm (26") Transparent (NC: 1304)	£82
DDU1202		609mm (24") Kuschall print (NC: 1304)	£82
DDU1204		559mm (22"), 609mm (24") Boa-Boa (NC: 1304)	£115
DDU1210		609mm (24") "Switzerland (NC: 1304)	£115
DDU1211		609mm (24") Pirates (NC: 1304)	£115
DDU1207		609mm (24") "Original" Carbon (NC: 1304)	£122
SAFETY			
DDU1300		Push brake, bent (NC: 1310)	◆
DDU1301		Push brake lever straight	○
DDU1302		Performance brake (NC: 1304, 1310)	£53
DDU1303		Scissor Brake (NC:1310) See brake table on back page for combinations with 609mm (24") rear wheels	£53
DDU1306		Pull brake	○


DRUM BRAKE		
DDU1304	Attendant drum brake	£259
DDU1310	Brake lever extension	£31
ANTI TIPPER Only 2 adaptions per side including Tubular arms, Anti tips, Tipper aid, Cane holder and Transit wheels		
DDU1702	Anti tipper right	£42
DDU1703	Anti tipper left	£42
DDU1704	Anti tipper both sides	£82
SEAT BELT		
DDU1730	Seat belt	£52
ACCESSORIES		
DDU1729	Rear frame protectors	£12
DDU1739	Front frame protectors	£23
DDU1745	Bottle holder and bottle	£11
DDU1746	Mobile phone holder	£12
DDU1708	Tipper-Aid, right (NC: 1706)	£30
DDU1709	Tipper-Aid, left (NC: 1707)	£30
CANE HOLDER		
DDU1706	Cane holder, right	£46
DDU1707	Cane holder, left	£46
TRANSIT WHEELS 3" NB: Maximum of 2 rear frame adaptions per side, including Tubular Armrests		
DDU1700	Transit wheels 3"	£147
MANUAL AIR PUMP		
DDU1715	Air pump	£9
TOOL KIT		
DDU1717	Tool kit	◆
PRICE		
Total Price		

Non Conforming Items Guide	
DDU0118	28cm Seat width
DDU0120	30cm Seat width
DDU0249	Rear wheel extension
DDU0267	5th position
DDU0268	6th position
DDU0283	Camber 3°
DDU0660	Flo-Tech Lite cushion
DDU0703	Height adjustable, rear set
DDU0712	Angle adjustable backrest, adjusted 102°
DDU0734	Fold down back (folding forward onto the seat)
DDU0750	Without push handles
DDU0820	6" Light wheel
DDU0821	6" Pneumatic Wheel
DDU0822	7" Light Wheel
DDU0904	24" Spider-Wheel
DDU0908	24" High Performance Wheel
DDU0910	Cross spoked
DDU0911	Radial spoked
DDU0913	24" Mountain Bike
DDU0915	22", 24", 26" Kuschall Starec
DDU0916	E-motion adaptation plate
DDU0917	24",25" Spinergy Wheel
DDU0918	25" High Performance Wheel
DDU1001	22", 24", 26" Light wheel 7 bar
DDU1002	22", 24" Profile tyre 7 bar
DDU1004	24" Solid tyre, profil, grey
DDU1006	24" Solid tyre (KIK type), black
DDU1007	22", 24" Solid, light wheel
DDU1020	Anti flat tyre tube protection
DDU1102	22", 24", 26" Supergripp, dist. 3 cm
DDU1104	24", 26" Steel, dist. 3 cm
DDU1106	24" Titanium handrim
DDU1120	24", 25" Surge LT handrim
DDU1300	Push brake bent
DDU1302	Performance brake
DDU1304	Attendant drum brake
DDU1310	Brake lever extension
DDU1405	Mud guards carbon, fixed
DDU1411	Tubular Armrest
DDU1416	Siderest height adjustable (with tools), removable, short
DDU1421	Siderest height adjustable (with tools), removable, long
DDU1430	Siderest flip-back, height adjustable, removable, short
DDU1431	Siderest flip-back, height adjustable, removable, long
DDU1440	Siderest with desk function, flip-back, height adjustable, removable, short
DDU1441	Siderest with desk function, flip-back, height adjustable, removable, long
DDU1491	Backrest stabilisation bar, swing away and removable
DDU1530	Heel loops
DDU1536	Insteepl loops
DDU1564	2 piece angle adjustable footplate set to front
DDU1565	2 piece angle adjustable footplate set back
DDU1702	Antitipper right
DDU1703	Antitipper left
DDU1706	Cane Holder Right
DDU1707	Cane Holder Left

**Adducted Frame,
Possible footwidth**
(SW - 14 cm)

- SB 32 cm : 18 cm
- SB 34 cm : 20 cm
- SB 36 cm : 22 cm
- SB 38 cm : 24 cm
- SB 40 cm : 26 cm
- SB 42 cm : 28 cm
- SB 44 cm : 30 cm
- SB 46 cm : 32 cm
- SB 48 cm : 34 cm
- SB 50 cm : 36 cm

Customer UL must be measured
taking away the cushion
thickness



**RSF / Rear
wheel size**

- 22" : 38-45 cm
- 24" : 39-48 cm
- 25" : 41-50 cm
- 26" : 42-50 cm

* these measurements can be achieved
as a special construction. Contact our
Specials department for further details

Folding Back Angle Possibilities		
Back Angle Choice	Adjustment Possibilities	
82°	86°	90°
86°	82°	90°
90°	82°	86°
94°	90°	98°
98°	90°	94°

Front wheel//Front wheel attachment	high mounted	standard mounted	low mounted
3"	38* - 47	41 - 50	44 - 53
4"	41 - 48	44 - 51	47 - 53
5"	44 - 49	47 - 52	50 - 53
6"	46 - 51	49 - 53	52 - 53
7"	48 - 52	51 - 53	54-55*

Scissor Brake compatibilities on a 24" Rear Wheel							
Scissor Brakes	Seat Depth	Rear Wheel Position					
		1	2	3	4	5	6
	32cm-38cm	N	N	N	N	Y	Y
	40cm	N	N	N	N	Y	Y
	42cm	N	N	N	Y	Y	Y
	44cm	N	N	Y	Y	Y	Y
	46cm	N	Y	Y	Y	Y	Y
	48cm	Y	Y	Y	Y	Y	Y
	50cm	Y	Y	Y	Y	Y	Y

Rear wheel extension
with an extra 2 positions



<Std axle Front

