



Yes, you can:

e-spare cat Now Available Online at
www.invacare.co.uk

Customer Name :	Account No :
Dealer :	Order No :
Address :	Order Qty :
Tel No :	Order Date :
Delivery address :	

Küschall® R33
DEALER
Maximum user weight 100kg (16 Stone)

Standard Features Active rigid wheelchair Highly configurable Standard light backrest upholstery Starec hub wheels Clothes guards, Calf strap 5 Different Seat Depths 6 Different Seat Widths Colour power set Standard Base Price £3,271
--

Please make your selection with a cross in the box ◆ Standard feature ○ Free of Charge

Customer weight =		Kg	
Code	X	Description	Dealer
DDR0080		Küschall® R33	£3,271
DDR0000		Special construction # Please fax or email the full requirements #	P.O.A.
FRAME ANGLE			
DDR0237		Küschall R33, frame angle 75° to the ground	◆
DDR0238		Küschall R33, frame angle 90° to the ground (NC:0811,0863,0813,0814,0818,0832,0817,0819)	○
SEATING			
DDR0626		Shock absorber (NC: 0925 with 0272)	◆
SHOCK ABSORBER USER WEIGHT			
DDR0635		Low Impact System - user weight < 70 kg	○
DDR0636		Low Impact System - user weight between 70 kg & 85 kg	○
DDR0637		Low Impact System - user weight between 85 kg & 100 kg	○
SEAT WIDTH (SW)			
DDR0124		Seat width 34 cm (NC: 0660)	○
DDR0126		Seat width 36 cm (NC: 0660)	○
DDR0128		Seat width 38 cm	○
DDR0130		Seat width 40 cm	○
DDR0132		Seat width 42 cm	○
DDR0134		Seat width 44 cm (NC:0660)	○
SEAT DEPTH (SD)			
DDR0414		Seat depth 37.5 cm (NC: 0660)	○
DDR0417		Seat depth 40 cm	◆
DDR0420		Seat depth 42.5 cm	○
DDR0423		Seat depth 45 cm (NC: 0660)	○
DDR0424		Seat depth 47.5 cm (NC: 0660)	○
SEAT POSITION (SP) Position 1 is more Active, Position 6 is less active			
DDR0263		154mm Position 1	○
DDR0264		131mm Position 2	○
DDR0265		108mm Position 3	○
DDR0266		85mm Position 4	◆
DDR0267		62mm Position 5	○
DDR0268		39mm Position 6	○
BACK HEIGHT (RH)			
DDR0548		RH 27 cm	○
DDR0580		RH 28,5 cm	○
DDR0550		RH 30 cm	○
DDR0503		RH 31,5 cm	○
DDR0506		RH 33 cm	○

Order No :

Customer Ref :

DDR0509		RH 34,5 cm		○
DDR0512		RH 36 cm		○
DDR0514		RH 37,5 cm		◆
DDR0516		RH 39 cm		○
DDR0518		RH 40,5 cm		○
DDR0520		RH 42 cm		○
DDR0522		RH 43,5 cm		○
DDR0524		RH 45 cm		○
DDR0526		RH 46,5 cm		○
FRONT SEAT-TO-FLOOR (SHv)				
DDR0848		FSF 46 cm	Available with 24" wheels only	○
DDR0849		FSF 47 cm	Available with 24" or 25" wheels	○
DDR0850		FSF 48 cm	Available with 24" or 25" wheels	○
DDR0851		FSF 49 cm	Available with 24" or 25" wheels	○
DDR0852		FSF 50 cm	Available with 24" or 25" wheels	◆
DDR0853		FSF 51 cm	Available with 25" wheels only	○
REAR SEAT-TO-FLOOR (SHh)				
Additional screws and spacers for 2 height adjustments are included with the chair				
DDR0948		RSF 43 cm	Available with 24" wheels only	◆
DDR0949		RSF 44 cm	Available with 24" or 25" wheels	○
DDR0950		RSF 45 cm	Available with 24" or 25" wheels	○
DDR0951		RSF 46 cm	Available with 24" or 25" wheels	◆
DDR0952		RSF 47 cm	Available with 24" or 25" wheels	○
DDR0953		RSF 48 cm	Available with 25" wheels only	○
KNEE-TO-HEEL LENGTH (UL) High Mount Footplate				
DDR1671		UL 30 cm		○
DDR1672		UL 31 cm		○
DDR1673		UL 32 cm		○
DDR1674		UL 33 cm		○
DDR1675		UL 34 cm		○
DDR1676		UL 35 cm		○
DDR1677		UL 36 cm		○
DDR1678		UL 37 cm		○
KNEE-TO-HEEL LENGTH (UL) Standard Footplate 75° frame: UL < SHv -3cm. 90° frame: UL < SHv -4cm				
DDR1680		UL 39 cm		○
DDR1681		UL 40 cm		○
DDR1682		UL 41 cm		○
DDR1683		UL 42 cm		○
DDR1684		UL 43 cm		◆
DDR1685		UL 44 cm		○
DDR1686		UL 45 cm		○
DDR1687		UL 46 cm		○
DDR1688		UL 47 cm		○
DDR1689		UL 48 cm		○
BACKREST				
BACKREST MATERIAL				
DDR0730		Aluminium		◆
DDR0731		Titanium	90g weight reduction	£112
BACK UPHOLSTERY				
DDR0722		Tension adjustable back, Black/grey		◆
DDR0721		Light upholstery	(Weight Saving 355gms)	£114
DDR0735		Aero upholstery with tension adjustable back		£158
DDR0746		Alcantara Special Edition 01	(Grey with Red Stripe)	£93
DDR0749		Alcantara Special Edition	(Black with Red Stripe)	£93
DDR0752		Alcantara Special Edition	(Grey with Black Stripe)	£93
DDR0748		Alcantara Black		£71
PUSH HANDLES				
DDR0750		Without push handles		◆
DDR0702		Standard push handles		£93
DDR0704		Foldable push handles		£141
SEATING				
SEAT UPHOLSTERY				
DDR0610		Standard		◆
DDR0615		Light seat	(Weight saving 60gms)	£39
DDR0611		Tension adjustable		£30

FOOTREST			
FOOTREST MATERIAL			
Titanium footrest reduces weight by 120g			
DDR1508		Aluminium	◆
DDR1509		Titanium	£77
FOOTREST (measure without cushion)			
Footrest flip to side not available up mounted			
DDR1510		Footrest tube, rigid	◆
DDR1570		Carbon footplate, angle adjustable	£141
SIDE GUARDS			
Clothes guards/Mud guards			
DDR1401		Carbon clothes guards	◆
DDR1403		Carbon mud guards (NC :1416,1421)	£141
ARMREST			
DDR1416		Siderest height adjustable, short	£68
DDR1421		Siderest height adjustable, long	£68
CUSHIONS AND CUSHION OPTIONS			
DDR0600		Seat cushion (NC: 0660)	£47
DDR0660		Flo Tech Lite cushion (Only for seat widths 38cm to 42 cm and seat depth 38cm to 44cm)	£101
DDR0603		Alcantara cushion Special Edition 01 (Grey with Red Stripe)	£65
DDR0607		Alcantara Special Edition (Black with Red Stripe)	£65
DDR0605		Alcantara cushion Black	£36
FRAME			
FRAME COLOUR, STANDARD			
DDR0300		Glossy black	○
DDR0314		Black Matt	○
DDR0354		Night shimmer	○
DDR0303		Royal red	○
DDR0302		Royal Blue	○
DDR0315		Silver grey	○
DDR0356		Ivory white	○
DDR0350		Purple fire	○
DDR0352		Emerald green	○
DDR0355		Midnight blue	○
DDR0357		Umber glaze	○
DDR0358		Smokey grey	○
DDR0359		Flashy yellow	○
DDR0363		Black-grey rustico	○
POLISHED FRAME			
DDR0310		Mirror natural	£188
COLOUR-POWER-SET Starec front castors included			
DDR0330		Colour-Power-Set techno-anthracite	○
DDR0331		Colour-Power-Set lemon	○
DDR0329		Colour-Power-Set comet red	○
DDR0336		Colour-Power-Set arctic blue	○
DDR0334		Colour-Power-Set b-blue	○
DDR0338		Colour-Power-Set chrome	○
FRONT WHEEL			
DDR0800		3" Supersport	○
DDR0803		4" Starec techno-anthracite	○
DDR0804		4" Starec lemon	○
DDR0808		4" Starec comet red	○
DDR0831		4" Starec arctic blue	○
DDR0807		4" Starec b-blue	○
DDR0809		4" Starec black	○
DDR0861		4" Starec chrome	○
DDR0863		5" Black rubber comfort wheel	◆
DDR0813		5" Starec techno anthracite	○
DDR0814		5" Starec lemon	○
DDR0818		5" Starec comet red	○
DDR0832		5" Starec arctic blue	○
DDR0817		5" Starec b-blue	○
DDR0819		5" Starec black	○
DDR0862		5" Starec chrome	○

CAMBER AXLE		
DDR0278	Carbon Axle	◆
CAMBER		
DDR0270	0° camber	○
DDR0271	3° camber	◆
DDR0272	6° camber (NC: 1302)	○
TETRA CLIP AXLE		
DDR0909	Tetra Clip	£60
REAR WHEEL SIZE		
DDR0902	24" wheel (NC: 0918)	○
DDR0925	25" wheel (NC: 0905,0908,0913)	○
REAR WHEEL TYPE		
DDR0905	24" Kuschall Starec wheel	◆
DDR0908	24" High performance, Black rim, Starec hub in CPS when selected (NC: 0910,0911,1004,1006,1104) No Spokeguards available	£177
DDR0918	25" High performance, Black rim, Starec hub in CPS when selected (NC:0910,0911,1001,1002,1004,1006,1102,1104,1106) No Spokeguards available	£258
DDR0913	24" Mountain Bike, fitted with specific 24" x 1.75" tyre, Aluminium or Supergrip handrim choice (NC: 0910,0911,1004,1006,1011,1010,1104,1106,1120,1202,1207,1209)	£159
DDR0917	24" Spinergy wheel (NC: 0910,0911,1302) No Spokeguards available for 24" or 25" Spinergys 25" Spinergy wheel (NC: 0910,0911,1302,1001,1002,1104,1006,1102,1104,1106)	£409
TYPE OF SPOKES		
DDR0910	Crossed spoked	◆
DDR0911	Radial spoked	£93
REAR WHEEL TYRE (No tyre choice for Mountain Bike Wheel)		
DDR1001	Light wheel 7 bar	○
DDR1002	Profile tyre 7 bar	◆
DDR1004	Solid tyre, profile, grey	£47
DDR1006	Solid tyre (KIK) black	£61
DDR1011	24", 25" Schwalbe speedrun	£60
DDR1010	Black Marathon Evolution tyre with kevlar lining 24" (NC: 1020)	£71
ANTI-FLAT		
DDR1020	Anti-flat tyre tube protection	£14
HANDRIM		
DDR1102	24" Supergrip, distance 3 cm	£114
DDR1104	24" Steel polished, distance 3 cm	£44
DDR1106	24" Titanium, distance 3 cm	£124
DDR1113	24", 25" Light Alu, distance 3 cm	○
DDR1114	609mm (24"), 635mm (25") Light Alu, distance 4 cm	◆
DDR1120	609mm (24"), 635mm (25") Surge LT 30mm dist	£230
SPOKE GUARDS		
DDR1200	24" transparent	£82
DDR1202	24" kusch	£82
DDR1207	24" "Original" Carbon	£115
DDR1209	24" R33 spoke guards	£82
SAFETY		
BRAKES		
DDR1300	Standard brake lever bent	◆
DDR1302	Performance brake	£53
DDR1303	Active brake Not possible on 25" rear wheels in combination with 6° camber & Low Impact System	£53
SAFETY OPTIONS		
ANTI TIPPER		
DDR1702	Anti-tipper tube right	£31
DDR1703	Anti-tipper tube left	£31
DDR1704	Anti-tipper tube both sides	£61
ACCESSORIES		
DDR1715	Air Pump	£9
DDR1739	Front frame protectors (Not compatible with upmounted footplate)	£23
DDR1717	Tool Kit	◆
DDR1745	Bottle holder with bottle	£11
DDR1746	Mobile phone holder	£12
PRICE		
Total price		£

Non Conforming Items Guide	
DDR0272	6° camber
DDR0660	Flo-tech cushion
DDR0811	5" Softroll
DDR0813	5" Starec techno anthracite
DDR0814	5" Starec lemon
DDR0817	5" Starec b-blue
DDR0818	5" Starec comet red
DDR0819	5" Starec black
DDR0832	5" Starec arctic blue
DDR0863	5" Black rubber comfort wheel
DDR0905	24" Kuschall starec wheel
DDR0908	24" high performance wheel
DDR0913	24" Mountain bike wheel
DDR0910	Crossed spoked
DDR0911	Radial spoked
DDR0918	25" High Performance wheel
DDR0925	25" wheel
DDR1001	24" Light wheel 7 bar
DDR1002	24" Profile tyre 7 bar
DDR1004	24" Solid tyre, profile, grey
DDR1006	24" Solid tyre black (KIK)
DDR1010	Marathon Evolution tyre
DDR1011	Schwalbe speedrun tyre
DDR1020	Anti-flat tyre tube protection
DDR1102	24" Supergrripp, distance 3 cm
DDR1104	24" Steel polished, distance 3 cm
DDR1106	24" Titanium, distance 3 cm
DDR1120	24" & 25" Surge Handrim
DDR1202	24" küsch spokeguard
DDR1207	24" "Original" Carbon spokeguard
DDR1209	24" Special R33 spoke guards
DDR1302	Performance Brake
DDR1416	Siderest height adjustable, short
DDR1421	Siderest height adjustable, long