

Invacare® Top End® XLT Handcycle - UK Price List and Order Form

For accuracy of ordering, please fax order to Invacare Customer Service on 01656 776220 or e-mail to: ordersuk@invacare.com



Provider name: _____ Account number: _____
 PO number: _____
 Phone number: _____ Fax number: _____
 Address: _____

Customer name: _____ Phone number: _____
 Address: _____

Ship to: Provider Customer (Order will be shipped to provider unless otherwise noted)

DATE:

QUOTE

ORDER

Please indicate your choice with a cross in this square


Indicate that this option has no upcharge

Manual Sport Wheelchair
 January 2012

Item	Description	Comments and /or non compatibilities information	Price
<input type="checkbox"/> XLT	Top End XLT		£1,899
FEATURES			
Standard	Aluminium frame / chrome-moly fork		<input type="radio"/>
Standard	Adjustable seat and back angle		<input type="radio"/>
Standard	Parking brake	Mounted on the right side of the handcycle	<input type="radio"/>
Standard	15° camber		<input type="radio"/>
Standard	Road crown compensator		<input type="radio"/>
SEATING OPTIONS			
Seat width			
<input type="checkbox"/>	Seat width 15" (app. 381 mm)		<input type="radio"/>
<input type="checkbox"/>	Seat width 17" (app. 432 mm)		<input type="radio"/>
BACK STYLE			
Standard	Narrow	Seat width minus 7", back height fixed at 20"	<input type="radio"/>
FOOTREST			
Standard	Adjustable footrest system	Multi-adjustable fore and aft with straps to securely cradle calves and feet. Leg guard attachments are recommended for persons who have no leg movement. The leg guard will protect the legs from getting into the chain.	<input type="radio"/>
CRANK SET			
Standard	Adjustable crank height		<input type="radio"/>
PROPULSION			
Standard	Shimano® 7 speed hub with reversing drum brake 175mm long x 12" wide crankset with one chain ring Ovalized ergonomic handpedals	7 speed front hub with internal hub brake and 7 speed twist shifter mounted on the left side of the frame but easily moveable to the right side. Adjustable crank height. <u>Braking system:</u> Shimano internal drum/coaster brake and Shimano HG50 Chain activated by pedalling backwards. <u>Gear shifter:</u> Shimano Tourney revoshifter indexed twist grip with gear display, 7 speeds. Check also quad options.	<input type="radio"/>
HAND PEDALS			
Standard	Ovalized ergonomic handpedals	Shape and angle is ergonomically designed for max power and are a lightweight ovalized aluminium shape with frictionless bearings.	<input type="radio"/>
WHEELS			
Standard	Quick release axles		<input type="radio"/>
Rear wheels type			
Standard	26" High Performance wheels	Rear hubs are precision heat treated, silver anodized with quick release stainless axles.	<input type="radio"/>
Tyres			
Standard	High pressure clincher tyres	26" x 1" High Performance clincher tyre / tube 100 psi	<input type="radio"/>
SAFETY AND COMFORT OPTIONS			
Standard	Seat and back cushions		<input type="radio"/>
Standard	Chain guard	Black plastic	<input type="radio"/>
Standard	Safety flag		<input type="radio"/>
Standard	Chest restraint: hook and loop		<input type="radio"/>
Standard	Reflectors		<input type="radio"/>
<input type="checkbox"/> XLT139	Quad twist shifter option	Recommended for customers with limited hand function, 3 pronged adaptation secured around the twist shifter	£50
<input type="checkbox"/> XLT22	Quad gloves	Quad gloves secure paraplegic's hand	£81
<input type="checkbox"/>	<input type="checkbox"/> Small		
<input type="checkbox"/>	<input type="checkbox"/> Medium		
<input type="checkbox"/>	<input type="checkbox"/> Large		
<input type="checkbox"/>	<input type="checkbox"/> X-Large		
<input type="checkbox"/> XLT8	Tri-pin quad handpedals	Tri-pin quad handles offer a 3 point contact and also allow customer to get on/off of handles independently, as to change gears, etc.	£237
<input type="checkbox"/> XLT84	Crank width adaptor set	Adds 1.5" per set, right and left specific.	£15
<input type="checkbox"/> XLT12	Mirror		£20
<input type="checkbox"/> XLT13	Water bottle & cage		£24
<input type="checkbox"/> XLT14	Cordless computer		£50

<input type="checkbox"/>	XLT15	Helmet		£50
<input type="checkbox"/>		Small (6 5/8" - 7")		
<input type="checkbox"/>		Medium (6 7/8" - 7 1/4")		
<input type="checkbox"/>		Large (7 1/4 - 7 5/8")		
<input type="checkbox"/>	XLT6	Front safety light		£17
<input type="checkbox"/>	XLT7	Rear flashing light		£17
<input type="checkbox"/>	XLT17	Crutch holder and strap		£84
<input type="checkbox"/>	XLT9	Leg guard attachment		£50
<input type="checkbox"/>	XLT19	Tow bar for wheelchair		£170
<input type="checkbox"/>	XLT28	Backpack hydration system		£66
<input type="checkbox"/>	TER126	Click-strap large	Pad length approx. 12"	£58
<input type="checkbox"/>	TER125	Click-strap medium	Pad length approx. 10"	£54
<input type="checkbox"/>	TER124	Click-strap small	Pad length approx. 8"	£51
<input type="checkbox"/>	RACE22	Top End t-shirt Size L		£9
<input type="checkbox"/>	RACE22	Top End t-shirt Size XL		£9
<input type="checkbox"/>	RACE13	Top End bike jersey Size L		£62
<input type="checkbox"/>	RACE13	Top End bike jersey Size XL		£62
UPHOLSTERY				
	Standard	Nylon upholstery with hook and loop tension adjustable seat sling. Seat restraint included.		○
FRAME COLOUR				
Contemporary colours				
<input type="checkbox"/>		Wet black		○
<input type="checkbox"/>		Black prism		○
<input type="checkbox"/>		White snow		○
<input type="checkbox"/>		Silver metallic		○
<input type="checkbox"/>		Stealth black		○
Way cool colours				
<input type="checkbox"/>		Electric blue		○
<input type="checkbox"/>		Midnight blue		○
<input type="checkbox"/>		Grasshopper green		○
Too hot colours				
<input type="checkbox"/>		Tangerine red		○
<input type="checkbox"/>		Electric red		○
<input type="checkbox"/>		Grape madness		○
<input type="checkbox"/>		Safety yellow		○
<input type="checkbox"/>		Sunkissed orange		○
<input type="checkbox"/>		Bubblegum pink		○

ADDITIONAL INFORMATION & MEASURING GUIDE



Material information
Frame: 2.5" x .0065 ovalized 6061 aluminium aircraft tubing
Frame design: single main tube, recumbent position
Fork assembly: 4130 Chrome-moly for strength stiffness. Raked angle for cornering stability.
Wheelbase: overall length 72"
Weight: overall weight 30-35 lbs (approx. 13.5 - 16 kg)
Turning radius: 18 feet for 360° circle, 12 feet with one backing manoeuvre
Static tip angle: 39° for 250 lb/113Kg subject
Transfer seat height: seat height front 13", seat height rear 9"
Back style / height: standard narrow back is ideal for competitive cyclists
Back angle: adjustable angle back
Steering: pivot steering with raked fork for controlled turning and stability. Road crown compensator is standard.
Warranty: 3 years on frame and fork to original owner/limited parts warranty
Shipping: shipped completely assembled except for rear wheels, shipping weight is 50 lbs (approx. 23 kg)

Invacare Ltd: www.invacare.co.uk

Pencoed Technology Park, Pencoed, Bridgend, CF35 5AQ T: 01656 776222 F: 01656 776220

All rights reserved. Trademarks are identified by the symbols ™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted.

Shimano and Rapid Fire are registered trademarks of Shimano Inc.

Delrin is registered trademark of E.I. duPont de Nemours and Company.

Specifications and prices are subject to change without notification.

©2011 Invacare Corporation - Form No. 96-270

