



January 2012

Yes, you can:

e-spare cat Now Available Online
at www.invacare.co.uk

Customer Name :	Account code:
Dealer :	Customer name:
Address :	Order no:
Tel No:	Order qty:
Delivery address :	Order date:

Küschall® Compact
DEALER
Maximum user weight 130kg (20.5 Stone)

Standard Features
<p>Lightweight aluminium folding frame 2 Camber Angles Clothes guards Starec Rear Wheels Aluminium handrims 2 Piece footplates 12 Different seat widths 10 Different seat depths 14 Frame Colours Standard Base Price £1,307</p>

Please make your selection with a cross in the box Standard feature Free of Charge

Code	X	Description	Dealer
DDP0012		Küschall® Compact	£1,307
DDP0000		Special Construction # Please fax or email the full requirements #	P.O.A.

SEATING CHOICE

FRONT FRAME

Important note: Push brake lever straight is not compatible with the standard frame option.

Code	Description	
DDP0211	Standard (NC: 1301)	◆
DDP0224	Abducted, Left (NC: 1300,1301)	○
DDP0205	Abducted, Right (NC: 1300,1301)	○
DDP0219	Abducted Left and Right (NC:1300)	○

SEAT WIDTH (SB) Note: Allow an extra 10mm on the actual seat width required.

Note: Seat widths 280mm - 320mm are only compatible with 1 or 2 piece carbon footplates

Code	Description	
DDP0118	Seat width 280 mm (NC:0660,1760,1118,1119) single cross brace Max user weight 130kg (20.5 Stone)	○
DDP0120	Seat width 300 mm (NC:0660,1760,1118,1119) single cross brace Max user weight 130kg (20.5 Stone)	○
DDP0122	Seat width 320 mm (NC:0660,1760,1118,1119) single cross brace Max user weight 130kg (20.5 Stone)	○
DDP0124	Seat width 340 mm (NC:0660,1760,1118,1119) single cross brace Max user weight 130kg (20.5 Stone)	○
DDP0126	Seat width 360 mm (NC:0660,1118,1119) single cross brace Max user weight 130kg (20.5 Stone)	○
DDP0128	Seat width 380 mm single cross brace Max user weight 130kg (20.5 Stone)	○
DDP0130	Seat width 400 mm single cross brace Max user weight 130kg (20.5 Stone)	○
DDP0132	Seat width 420 mm single cross brace Max user weight 130kg (20.5 Stone)	○
DDP0134	Seat width 440 mm (NC: 0660) double cross brace Max user weight 130kg (20.5 Stone)	○
DDP0136	Seat width 460 mm (NC: 0660) double cross brace Max user weight 130kg (20.5 Stone)	○
DDP0138	Seat width 480 mm (NC: 0660) double cross brace Max user weight 130kg (20.5 Stone)	○
DDP0140	Seat width 500 mm (NC: 0660) double cross brace Max user weight 130kg (20.5 Stone)	○

SEAT DEPTH (ST)

Code	Description	
DDP0408	Seat depth 320 mm (NC: 0660)	○
DDP0410	Seat depth 340 mm (NC: 0660)	○
DDP0412	Seat depth 360 mm (NC: 0660)	○
DDP0415	Seat depth 380 mm	○
DDP0417	Seat depth 400 mm	◆
DDP0420	Seat depth 420 mm	○
DDP0422	Seat depth 440 mm	○
DDP0424	Seat depth 460 mm (NC: 0660)	○
DDP0427	Seat depth 480 mm (NC: 0660)	○
DDP0429	Seat depth 500 mm (NC: 0660)	○

BACK HEIGHT (RH)			
DDP0550		RH 300 mm (Straight fixed backposts supplied as standard) (NC: 0734,1491)	○
DDP0503		RH 315 mm (Straight fixed backposts supplied as standard) (NC: 0734,1491)	○
DDP0506		RH 330 mm (Straight fixed backposts supplied as standard) (NC: 0734,1491)	○
DDP0509		RH 345 mm (NC: 1491)	○
DDP0512		RH 360 mm (NC: 1491)	○
DDP0514		RH 375 mm (NC: 1491)	○
DDP0516		RH 390 mm (NC: 1491)	○
DDP0518		RH 405 mm	○
DDP0520		RH 420 mm	◆
DDP0522		RH 435 mm	○
DDP0524		RH 450 mm	○
DDP0526		RH 465 mm	○
DDP0528		RH 480 mm	○
DDP0530		RH 495 mm	○
DDP0532		RH 510 mm	○
FRONT SEAT-TO-FLOOR HEIGHT (SHv)			
DDP0840		FSF 380 mm	○
DDP0841		FSF 390 mm	○
DDP0842		FSF 400 mm	○
DDP0843		FSF 410 mm	○
DDP0844		FSF 420 mm	○
DDP0845		FSF 430 mm	○
DDP0846		FSF 440 mm	○
DDP0847		FSF 450 mm	○
DDP0848		FSF 460 mm	○
DDP0849		FSF 470 mm	○
DDP0850		FSF 480 mm	○
DDP0851		FSF 490 mm	○
DDP0852		FSF 500 mm	◆
DDP0853		FSF 510 mm	○
DDP0854		FSF 520 mm	○
DDP0855		FSF 530 mm	○
REAR SEAT-TO-FLOOR HEIGHT (SHh)			
Mudguards in combination with a rear wheel extension are not compatible with rear seat to floor heights of 410mm or lower			
DDP0943		RSR 380 mm (NC: 1405 + 0249)	○
DDP0944		RSR 390 mm (NC: 1405 + 0249)	○
DDP0945		RSR 400 mm (NC: 1405 + 0249)	○
DDP0946		RSR 410 mm (NC: 1405 + 0249)	○
DDP0947		RSR 420 mm	○
DDP0948		RSR 430 mm	○
DDP0949		RSR 440 mm	○
DDP0950		RSR 450 mm	◆
DDP0951		RSR 460 mm	○
DDP0952		RSR 470 mm	○
DDP0953		RSR 480 mm	○
DDP0954		RSR 490 mm	○
DDP0955		RSR 500 mm	○
KNEE-TO-HEEL LENGTH (UL) Please note that on knee to heel lengths less than 350mm and depending configuration, the footplate may be upmounted! On 70° frame, knee to heel = (SHv - 20mm), On 80° & 90° frames, knee to heel = (SHv - 30mm)			
DDP1662		UL 210 mm (NC: 1562,1563,1564,1565)	○
DDP1663		UL 220 mm (NC: 1562,1563,1564,1565)	○
DDP1664		UL 230 mm (NC: 1562,1563,1564,1565)	○
DDP1665		UL 240 mm (NC: 1564,1565) (NC: 1540/1542 & 1562/1563)	○
DDP1666		UL 250 mm (NC: 1564,1565) (NC: 1540/1542 & 1562/1563)	○
DDP1667		UL 260 mm (NC: 1564,1565) (NC: 1540/1542 & 1562/1563)	○
DDP1668		UL 270 mm (NC: 1564,1565)	○
DDP1669		UL 280 mm (NC: 1564,1565)	○
DDP1670		UL 290 mm (NC: 1564,1565)	○
DDP1671		UL 300 mm (NC: 1564,1565)	○
DDP1672		UL 310 mm (NC: 1564,1565)	○
DDP1673		UL 320 mm (NC: 1564,1565)	○
DDP1674		UL 330 mm (NC: 1540 & 1564/1565)	○
DDP1675		UL 340 mm	○
DDP1676		UL 350 mm	○
DDP1677		UL 360 mm	○
DDP1678		UL 370 mm	○
DDP1679		UL 380 mm	○
DDP1680		UL 390 mm	○
DDP1681		UL 400 mm	○

DDP1682		UL 410 mm	○
DDP1683		UL 420 mm	○
DDP1684		UL 430 mm	◆
DDP1685		UL 440 mm	○
DDP1686		UL 450 mm	○
DDP1687		UL 460 mm	○
DDP1688		UL 470 mm	○
DDP1689		UL 480 mm	○
DDP1690		UL 490 mm	○
DDP1691		UL 500 mm	○
DDP1692		UL 510 mm	○
BACKREST CHOICE			
BACK UPHOLSTERY			
DDP0700		Standard, light upholstery (NC: 0750,0703,1491,1731,1732,1733)	◆
DDP0722		Tension velcro adjustable back black/grey with cover to the seat	£97
DDP0735		Aero back, tension velcro adjustable	£158
PUSH HANDLES			
DDP0702		Push-handles fixed, long	◆
DDP0750		without push handles	○
DDP0751		Push handles fixed, mini (NC: 1304)	○
DDP0703		Height adjustable, rear set	£180
DDP0733		Height adjustable, integrated	£87
DDP0704		Foldable push handles. (NC: 1304,1491)	£140
BACK TUBES <i>If no option is selected, straight backrest tubes will be assembled</i>			
DDP0705		Lumbal tubes (NC: 0734)	○
DDP0720		Angle adjustable back (NC: 0734) (Backrest angle must be selected from the 6 options below)	£143
DDP0734		Fold down back (folding forward onto the seat) (NC: 0712,1411) (Backrest angle must be selected, see table on back page for angle possibilities)	£52
BACKREST ANGLE			
DDP0707		Angle adjustable backrest, adjusted 82°	○
DDP0718		Angle adjustable backrest, adjusted 86°	○
DDP0709		Angle adjustable backrest, adjusted 90°	○
DDP0719		Angle adjustable backrest, adjusted 94°	○
DDP0711		Angle adjustable backrest, adjusted 98°	○
DDP0712		Angle adjustable backrest, adjusted 102°	○
BACKBRACE			
DDP1491		Backbrace swing away and removable	£88
SEATING			
DDP0611		Standard adjustable seat upholstery with 1 central velcro band	◆
LEGRESTS Note: Legrests from the Action range are compatible with seat widths of 360mm-500mm			
DDP1540		Legrests 70°	◆
DDP1541		Legrests 80° (NC: 0822)	○
DDP1542		Legrests 90° (NC: 0822)	○
FOOTREST For Seat Width 280-320mm only 1pc or 2pc Carbon footplates are possible			
DDP1562		2 pc. plastic, set to front	◆
DDP1563		2 pc. plastic, set to back	○
DDP1564		2 pc. plastic angle adjustable set to front	£92
DDP1565		2 pc. plastic angle adjustable set back	£92
DDP1578		2 pc. carbon footplates, set to the front (NC: 1530,1536)	£237
DDP1579		2 pc. carbon footplates, set to the rear (NC: 1530,1536)	£237
DDP1566		1 pc. carbon Flip-up, adjustable, fixed on 1 side when folded, set to front (NC: 1530,1536)	£118
DDP1567		1 pc. carbon Flip-up, adjustable, fixed on 1 side when folded, set to back (NC: 1530,1536)	£118
FOOT AND LEG SUPPORT			
DDP1534		Calf strap	○
DDP1530		Heel loops	○
DDP1536		Instep loops	○
SIDE GUARDS			
CLOTHES & MUD GUARDS			
DDP1400		Clothes guards carbon, fixed (NC:1416,1421,1430,1431,1440,1441)	◆
DDP1405		Mud guards carbon, fixed (NC: 1416,1421,1430,1431,1440,1441)	£112
DDP1408		Mud guards carbon, removable (NC:1416,1421,1430,1431,1440,1441) (Adds 20mm to the total width of the chair)	£199

ARMRESTS (Includes clothesguards) (Adds 20mm to total width of the chair)		
DDP1416	Armrest height adjustable (with tools), removable, short	○
DDP1421	Armrest height adjustable (with tools), removable, long	○
DDP1430	Armrest flip-back height adjustable removable, short (in frame colour, with carbon clothesguards) (NC: 1731,1732,1733,1760)	£293
DDP1431	Armrest flip-back height adjustable removable, long (in frame colour, with carbon clothesguards) (NC: 1731,1732,1733,1760)	£293
DDP1440	Armrest with desk function, flip-back, height adjustable, removable, short (NC: 1731,1732,1733,1300,1302,1760,0283)	£150
DDP1441	Armrest with desk function, flip-back, height adjustable, removable, long (NC: 1731,1732,1733,1300,1302,1760,0283)	£180
DDP1411	Tubular armrests height adjustable,removable,swingaway, comfort padded armpad, no clothesguards	£170
If Tubular armrests are selected, only 1 other adaption can be fitted. Adaptions include Anti tips, Tipper tube, Cane holder and Transit wheels		
SEATING OPTIONS		
SEAT CUSHION		
DDP0600	Standard seat cushion	£47
DDP0660	Flo Tech Lite seat cushion. (Only available for seat width 380 to 420mm and seat depth 380 to 440mm)	£101
DDP0606	Seat Bag, (click in and removable)	£36
COLOURS		
FRAME		
DDP0300	Glossy-Black	○
DDP0314	Black Matt	○
DDP0356	Ivory White	○
DDP0302	Royal Blue	○
DDP0303	Royal Red	○
DDP0350	Purple Fire	○
DDP0315	Silver Grey	○
DDP0354	Night Shimmer	○
DDP0352	Emerald Green	○
DDP0355	Midnight Blue	○
DDP0357	Umber Glaze	○
DDP0358	Smokey Grey	○
DDP0359	Flashy Yellow	○
DDP0363	Black-grey rustico	○
COLOUR-POWER-SET		
DDP0330	Colour-Power-Set techno-anthracite	£158
DDP0331	Colour-Power-Set lemon	£158
DDP0332	Colour-Power-Set apricot	£158
DDP0334	Colour-Power-Set b-blue	£158
DDP0329	Colour-Power-Set comet red	£158
DDP0337	Colour-Power-Set copper glaze	£158
DDP0336	Colour-Power-Set artic blue	£158
DDP0338	Colour-Power-Set chrome	£158
FRONT WHEEL ATTACHMENT		
DDP0864	High mounted (Frame will be 30mm lower)	○
DDP0865	Standard mounted	◆
DDP0866	Low mounted (Frame will be 30mm higher)	○
FRONT WHEEL		
DDP0800	75mm (3") Supersport	○
DDP0801	100mm (4") Everyday	£45
DDP0803	100mm (4") Starec techno-anthracite	£45
DDP0804	100mm (4") Starec lemon	£45
DDP0805	100mm (4") Starec apricot	£45
DDP0807	100mm (4") Starec b-blue	£45
DDP0808	100mm (4") Starec comet red	£45
DDP0830	100mm (4") Starec copper glaze	£45
DDP0831	100mm (4") Starec artic blue	£45
DDP0809	100mm (4") Starec black	£45
DDP0861	100mm (4") Starec chrome	£45
DDP0824	100mm (4") Skater Wheel	£45
DDP0810	125mm (5") Light Wheel	○
DDF0863	125mm (5") Comfort Wheel black rubber	◆
DDP0813	125mm (5") Starec techno anthracite	£45
DDP0814	125mm (5") Starec lemon	£45
DDP0815	125mm (5") Starec apricot	£45
DDP0817	125mm (5") Starec b-blue	£45
DDP0818	125mm (5") Starec comet red	£45
DDP0833	125mm (5") Starec copper glaze	£45
DDP0832	125mm (5") Starec artic blue	£45
DDP0819	125mm (5") Starec black	£45
DDP0862	125mm (5") Starec chrome	£45
DDP0820	150mm (6") Light wheel	○
DDP0821	150mm (6") Pneumatic Wheel	£18
DDP0822	175mm (7") Light Wheel	○

REAR WHEEL CHOICE		
ONE-ARM-DRIVE SYSTEM (Complete system with Pneumatic profile grey tyres & aluminium handrims)		
DDP1118	One-arm-drive Right side (NC: 1304,0909,0283,0266,0267,0268) (Not possible with spokeguards)	£617
DDP1119	One-arm-drive Left side (NC: 1304,0909,0283,0266,0267,0268) (Not possible with spokeguards)	£617
REAR WHEEL ADAPTION PLATE		
DDP0250	Standard fixation (NC: 0267,0268)	◆
DDP0249	Rear wheel extension (Black in colour)	£47
REAR WHEEL POSITION		
DDP0263	1st position (NC: 1304)	○
DDP0264	2nd position	○
DDP0265	3rd position	◆
DDP0266	4th position	○
DDP0267	5th position	○
DDP0268	6th position	○
QUAD QUICK RELEASE AXLE		
DDP0909	Tetra Quick release axle (NC: 1304)	£60
REAR WHEEL SIZE		
DDP0901	559mm (22") (NC: 0904,0908,0918,0913,0917,0916)	○
DDP0902	609mm (24")	◆
DDP0925	635mm (25") (NC: 0916,0904,0915,0908,0913,0283,1102,1104)	○
DDP0903	660mm (26") (NC: 0283,0916,0904,0908,0918,0913,0917)	○
REAR WHEEL CAMBER (3° Not available with drum brakes and performance brakes combined)		
DDP0282	Camber 0°	◆
DDP0283	Camber 3° (NC: 1302,1304, 0916)	○
REAR WHEEL TYPE		
DDP0904	609mm (24") Spider-Wheel (NC: 0910,0911,1104,1120) No spokeguards available	£158
DDP0915	559mm (22"), 609mm (24"), 660mm (26") Kuschall Starec (Hub in CPS if selected)	◆
DDP0908	609mm (24") High Performance (Starec hub) (NC: 0910,0911,1004,1006,1007,1104) No spokeguards available	£177
DDP0918	635mm (25") High Perf, with Marathon tyre & Alum handrim (Starec hub) No spokeguards available	£258
DDP0913	609mm (24") Mountain Bike with 24" x 1.75" specific tyre. (NC:0910,0911,1104,1106,1120) Available with Transparent or Boa-Boa spokeguards only	£159
DDP0917	609mm (24") Spinergy wheel (NC: 0910,0911,1104) No spokeguards available 635mm (25") Spinergy wheel (NC: 0910,0911,1001,1002,1004,1006,1007,1102,1104,1106)	£409
TYPES OF SPOKES (Select only with Kuschall Starec wheels)		
DDP0910	Cross spoked	○
DDP0911	Radial spoked	○
E-MOTION ADAPTATION (Does not include e-motion wheels) Only compatible with 609mm (24") wheels & 0 Degree camber		
DDP0916	E-motion adaptation plate (NC: 1702,1703,1706,1707)	£182
REAR WHEEL TYRE		
DDP1001	559mm (22"), 609mm (24"), 660mm (26") Light wheel 7 bar	◆
DDP1002	559mm (22"), 609mm (24") Profil tyre 7 bar	○
DDP1004	609mm (24") Solid tyre, profil, grey	£47
DDP1006	609mm (24") Solid tyre (KIK type), black	£61
DDP1007	559mm (22"), 609mm (24") Solid, light wheel	£47
DDP1010	609mm (24"), 635mm (25"), 660mm (26") Marathon Evolution tyre from Schwalbe (NC: 1020)	£71
DDP1011	609mm (24"), 635mm (25") Schwalbe speedrun	£60
ANTI-FLAT		
DDP1020	Anti Flat tyre tube protection	£14
HANDRIM		
DDP1113	559mm (22"), 609mm (24"), 635mm (25"), 660mm (26") Light Alu, dist. 3 cm	◆
DDP1114	559mm (22"), 609mm (24"), 635mm (25"), 660mm (26") Light Alu, dist. 4 cm	○
DDP1102	559mm (22"), 609mm (24"), 660mm (26") Supergripp, dist. 3 cm	£114
DDP1104	609mm (24"), 660mm (26") Steel, dist. 3 cm	£44
DDP1106	609mm (24") Titanium, dist. 3 cm	£124
DDP1120	609mm (24"), 635mm (25") Surge LT 30mm handrim, 3cm distance	£230
SPOKE GUARDS		
DDP1200	559mm (22"), 609mm (24"), 660mm (26") Transparent (NC: 1304)	£82
DDP1202	609mm (24") Kuschall print (NC: 1304)	£82
DDP1204	559mm (22"), 609mm (24") Boa-Boa (NC: 1304)	£115
DDP1207	609mm (24") "Original" Carbon (NC: 1304)	£122
DDP1210	609mm (24") "Switzerland (NC: 1304)	£115
DDP1211	609mm (24") Pirates (NC: 1304)	£115

SAFETY			
BRAKES			
DDP1306		Pull brake	◆
DDP1300		Push brake, bent (NC: 1310)	○
DDP1301		Push brake lever straight	○
DDP1302		Performance Brake (NC :1304, 1310)	£53
DDP1303		Scissor Brake (NC:1310) See brake table on back page for combinations with 609mm (24") rear wheels	£53
DRUM BRAKE			
DDP1304		Attendant drum brake	£259
BRAKE LEVER EXTENSION			
DDP1310		Brake extension lever	£31
ANTI TIPPER <i>Only 2 adaptions per side including Tubular arms, Anti tips, Tipper aid, Cane holder and Transit wheels</i>			
DDP1702		Anti tipper right	£42
DDP1703		Anti tipper left	£42
DDP1704		Anti tipper both sides	£82
SEAT BELT			
DDP1730		Seat belt	£52
ACCESSORIES			
DDP1729		Rear frame protector	£12
DDP1731		Lateral Support, both sides	£315
DDP1732		Lateral Support, right	£158
DDP1733		Lateral Support, left	£158
DDP1760		Therapy Table	£177
DDP1745		Bottle holder and bottle	£11
DDP1746		Mobile phone holder	£12
TIPPER AID			
DDP1708		Tipper-Aid, right (NC: 1706)	£30
DDP1709		Tipper-Aid, left (NC: 1707)	£30
CANE HOLDER			
DDP1706		Cane holder, right	£46
DDP1707		Cane holder, left	£46
TRANSIT WHEELS 75mm (3")			
DDP1700		Transit wheels 75mm (3")	£147
MANUAL AIR PUMP			
DDP1715		Air pump	£9
TOOL KIT			
DDP1717		Tool kit	◆
PRICE			
Total Price			£

Non Conforming Items Guide	
DDP0249	Rear wheel extension
DDP0266	4th position
DDP0267	5th position
DDP0268	6th position
DDP0283	Camber 3°
DDP0660	Flo-Tech Lite seat cushion
DDP0703	Height adjustable rear set
DDP0712	Angle adjustable backrest, adjusted 102°
DDP0734	Foldable back (to the front)
DDP0750	without push handles
DDP0822	7" Light Wheel
DDP0904	24" Spider-Wheel
DDP0908	24" High Performance (Starec hub)
DDP0909	Quad quick release axles
DDP0910	Cross spoked
DDP0911	Radial spoked
DDP0913	24" Mountain Bike
DDP0915	22", 24", 26" Kuschall Starec (Hub in CPS if selected)
DDP0916	E-motion adaptation plate
DDP0917	24", 25" Spinergy wheel
DDP0918	25" High Performance (Starec hub)
DDP1001	22", 24", 26" Light wheel 7 bar
DDP1002	22", 24" Profile tyre 7 bar
DDP1004	24" Solid tyre, profil, grey
DDP1006	24" Solid tyre (KIK type), black
DDP1007	22", 24" Solid, light wheel
DDP1020	Anti-flat tyre tube protection
DDP1102	22", 24", 26" Supergripp, dist. 3 cm
DDP1104	24", 26" Steel, dist. 3 cm
DDP1106	24" Titanium, dist. 3 cm
DDP1118	One arm drive system RH
DDP1119	One arm drive system LH
DDP1120	Surge LT handrim
DDP1300	Push brake bent
DDP1301	Push brake lever straight
DDP1302	Performance brake
DDP1304	Attendant drum brake
DDP1310	Brake extension lever
DDP1405	Mud Guards Carbon fixed
DDP1411	Tubular Armrests
DDP1416	Siderest height adjustable (with tools), removable, short
DDP1421	Siderest height adjustable (with tools), removable, long
DDP1430	Siderest flip-back, height adjustable, removable, short
DDP1431	Siderest flip-back, height adjustable, removable, long
DDP1440	Siderest with desk function, flip-back, height adjustable, removable, short
DDP1441	Siderest with desk function, flip-back, height adjustable, removable, long
DDP1491	Backbrace swingaway and removable
DDP1530	Heel loops
DDP1536	Insteep loops
DDP1540	Legrest, angle open, 70°
DDP1542	Legrest, angle closed, 90°
DDP1562	2 pc. plastic, set to the front
DDP1563	2 pc. plastic, set back
DDP1564	2 pc. angle adjustable set to the front
DDP1565	2 pc. angle adjustable, set back
DDP1702	Anti tipper right
DDP1703	Anti tipper left
DDP1706	Cane holder, right
DDP1707	Cane holder, left
DDP1731	Lateral Support, both sides
DDP1732	Lateral Support, right
DDP1733	Lateral Support, left
DDP1760	Therapy Table adjustable

Folding Back Angle Possibilities		
Back Angle Choice	Adjustment Possibilities	
82°	86°	90°
86°	82°	90°
90°	82°	86°
94°	90°	98°
98°	90°	94°

* these measurements can be achieved as a special construction. Contact our Specials department for further details

Front wheel/Front wheel attachment	high mounted	standard mounted	low mounted
3"	38* - 47	41 - 50	44 - 53
4"	41 - 48	44 - 51	47 - 53
5"	44 - 49	47 - 52	50 - 53
6"	46 - 51	49 - 53	52 - 53
7"	48 - 52	51 - 53	54-55*




Rear wheel extension offers positions 5 and 6

Scissor Brake compatibilities on a 24" Rear Wheel							
Scissor Brakes	Seat Depth	Rear Wheel Position					
		1	2	3	4	5	6
	32cm-38cm	N	N	N	N	Y	Y
	40cm	N	N	N	N	Y	Y
	42cm	N	N	N	Y	Y	Y
	44cm	N	N	Y	Y	Y	Y
	46cm	N	Y	Y	Y	Y	Y
	48cm	Y	Y	Y	Y	Y	Y
	50cm	Y	Y	Y	Y	Y	Y

CPS-Set Includes :

Starec-hub, forks, adapter plates, 22" and 24" Starec rear wheel rims, Quick release axle caps for Starec hubs, bearing block, seatprofile caps

Customer UL must be measured taking away the cushion thickness



SHh / Rear wheel size

22" : 38-45 cm

24" : 39-48 cm

25" : 41-50 cm

26" : 42-50 cm

<Std axle Front

