

Invacare Top End Force R Handcycle - UK Price List and Order Form

For accuracy of ordering, please fax order to Invacare Customer Service on 01656 776220 or e-mail to: ordersuk@invacare.com

www.topendwheelchair.co.uk



Provider name: _____ Account number: _____
 PO number: _____
 Phone number: _____ Fax number: _____
 Address: _____

Customer name: _____ Phone number: _____
 Address: _____

Ship to: Provider Customer (Order will be shipped to provider unless otherwise noted)



DATE:

QUOTE

ORDER

Manual Sport Wheelchair
 January 2012


- Please indicate your choice with a cross in this square
 Indicates that this option has no upcharge

Item	Description	Comments and/or non compatibilities information	Price
<input type="checkbox"/> FRCR	Top End Force R Handcycle		£4,424
<input type="checkbox"/> FRCR5	Force R Handcycle less front and rear wheels		£4,270
SEATING OPTIONS			
Seat width			
<input type="checkbox"/>	Seat width 13"	Inside lower frame 	<input type="radio"/>
<input type="checkbox"/>	Seat width 14"		<input type="radio"/>
<input type="checkbox"/>	Seat width 15"		<input type="radio"/>
<input type="checkbox"/>	Seat width 16"		<input type="radio"/>
<input type="checkbox"/>	Seat width 17"		<input type="radio"/>
<input type="checkbox"/>	Seat width 18"		<input type="radio"/>
Wheel base			
<input type="checkbox"/>	Standard	Wheelbase 60"	<input type="radio"/>
<input type="checkbox"/>		Wheelbase 62"	<input type="radio"/>
FOOTREST OPTIONS			
Multi adjustable fore and aft with straps to securely cradle calves and feet. Straps are 4" wide and feature safety straps for ankles.			
Lower fork style			
<input type="checkbox"/>	Standard	Aerodynamic narrow	Non compatible with disc brakes
<input type="checkbox"/>		Wide style	Required for disc brakes
Footrest width			
<input type="checkbox"/>		3 3/4" (approx. 95 mm)	
<input type="checkbox"/>		4 3/4" (approx. 120 mm)	
<input type="checkbox"/>	XLT83	Connected footrest with bumper	Not recommended for users under 5'8" tall - Can be custom built upon request for shorter leg lengths
			£55
COMPONENTS CHOICES: SHIFTER / BRAKE / FRONT WHEEL DERAILLEUR / CASSETTES PACKAGES			
Choose either 9 speed standard or 10 speed option and then indicate all the information respectively for component choice.			
<input type="checkbox"/> Standard: 9 speed			
<input type="checkbox"/>		SRAM® 9 speed Grip Shifter and SRAM® X.O. Derailleur / SRAM® Hollow Pin Chain	<input type="radio"/>
	Mounting position: Shifter / Brake Lever		
<input type="checkbox"/>		Right handle (standard)	<input type="radio"/>
<input type="checkbox"/>		Left handle	<input type="radio"/>
<input type="checkbox"/>	XLT64	SRAM® 9 speed Trigger Shifter and SRAM® X.O. Derailleur / SRAM® Hollow Pin Chain	Thumb and forefinger operated Super light weight, high quality
	Mounting position: Shifter		
<input type="checkbox"/>		Right handle mount	Left handle mount is not possible
<input type="checkbox"/>		Right fork mount	Left handle mount is not possible
	Mounting position: Brake Lever		
<input type="checkbox"/>		Right handle mount	<input type="radio"/>
<input type="checkbox"/>		Right fork mount	<input type="radio"/>
<input type="checkbox"/>		Left handle mount	<input type="radio"/>
<input type="checkbox"/>		Left fork mount	<input type="radio"/>
Front wheel cassette choice for 9 speed			
<input type="checkbox"/>	Std	SRAM® 11-34 (standard)	<input type="radio"/>
<input type="checkbox"/>	XLT58	SRAM® 11-23	Cassette with closer gear ratios - Not best for hilly terrain
<input type="checkbox"/>	XLT59	SRAM® 11-26	Cassette with closer gear ratios - Not best for hilly terrain
			£33
			£33

<input type="checkbox"/> Option: 10 speed				
<input type="checkbox"/>	XLT79	SRAM® XX 10 speed Shifter and brake SRAM® X.X. Derailleur / SRAM® Hollow Pin Chain	Right handle mounted shifter and brake Fork mounted locking second brake	£598
<input type="checkbox"/>	XLT80	SRAM® Red 10 speed Shifter and brake SRAM® X.X. Derailleur / SRAM® Hollow Pin Chain	Fork mounted shifter and brake Choose position below for second brake	£598
Only for XLT80: Second brake mounting position				
<input type="checkbox"/>		Left handle non-locking (standard)		<input type="radio"/>
<input type="checkbox"/>		Derailleur post non-locking	Only for manual derailleur	<input type="radio"/>
<input type="checkbox"/>		Fork mounted locking		<input type="radio"/>
Front wheel cassette choice for 10 speed				
<input type="checkbox"/>	Std	SRAM® 11-36 (standard)		<input type="radio"/>
<input type="checkbox"/>	XLT82	SRAM® 11-28	Cassette with closer gear ratios - Not best for hilly terrain	£33
FORK				
<input type="checkbox"/>	Standard	Fixed length / Vertical crank height adjustment	Rigid and lightest choice	<input type="radio"/>
<input type="checkbox"/>	XLT81	Fixed length and fixed crank height	Drawing required -see page 4	£222
<input type="checkbox"/>	XLT49	Diagonal and vertical crank height adjustment	Heavier than standard, less rigid	£201
HANDPEDAL SIZE				
Handpedals are made of lightweight aluminium and are ovalized to precisely fit the hand. The bearing which allow the pedal to rotate while pedaling are located in the handle and not in the crank arm, which makes it very comfortable.				
<input type="checkbox"/>	Standard	Medium handpedal size (1 ¾")		<input type="radio"/>
<input type="checkbox"/>	XLT72	Small handpedal size (1 ½")		£50
<input type="checkbox"/>	XLT73	Large handpedal size (2 ¼")		£50
<input type="checkbox"/>	XLT56	Power plate upgrade	Flat aluminium plates welded at the lower end of the pedal	£33
CRANK ARMS: TAPERED CRANK (choose length and width)				
Top End V crank set with integral pedal bearings are standard. Length of crank set depends on the length of the arms. To choose the width for the crank, consider the end user's shoulder width.				
Length in mm				
<input type="checkbox"/>		155 mm	Length is measured perpendicular from the centre of the crank to the centre of the hand pedal spindle.	<input type="radio"/>
<input type="checkbox"/>		160 mm		<input type="radio"/>
<input type="checkbox"/>		165 mm		<input type="radio"/>
<input type="checkbox"/>	Standard	170 mm		<input type="radio"/>
<input type="checkbox"/>		175 mm		<input type="radio"/>
<input type="checkbox"/>		180 mm		<input type="radio"/>
Width in inches				
<input type="checkbox"/>		12" (305 mm)	The width is measured horizontally centre to centre to the point at the middle of the hand pedal.	<input type="radio"/>
<input type="checkbox"/>		13" (330 mm)		<input type="radio"/>
<input type="checkbox"/>	Standard	14" (356 mm)		<input type="radio"/>
<input type="checkbox"/>		15" (381 mm)		<input type="radio"/>
<input type="checkbox"/>		16" (406 mm)		<input type="radio"/>
Extra crank arms				
<input type="checkbox"/>	XLT46	Length: _____ Width: _____		£99
CHAIN RING CHOICE (operated by crank mounted shifter)				
<input type="checkbox"/>	Standard	48-39-24		<input type="radio"/>
<input type="checkbox"/>	XLT41	Chain ring option 53-39-28		£66
<input type="checkbox"/>	XLT47	Rotor® Q-Ring option 53-39-28	28 chain will not be Q-ring	£380
PERFORMANCE UPGRADE				
<input type="checkbox"/>	XLT60	SRM® Training System - Professional road State of the art measuring tool	Not available with 3 chainrings or a chainring less than 39 teeth. Power meter is integrated directly into the cranks.	£2,575
BRAKE CHOICE				
<input type="checkbox"/>	Standard	Caliper-Cane Creek® SCR-3		<input type="radio"/>
<input type="checkbox"/>	XLT55	Caliper-Cane Creek® SL with titanium hardware	Ultra light, with titanium nuts and bolts	£102
<input type="checkbox"/>	XLT66	Disc brake	Only available with wide style fork and Std or Corima wheels	£187
WHEELS				
	Standard	Threaded axles only		<input type="radio"/>
Camber				
<input type="checkbox"/>		Camber 0°		<input type="radio"/>
<input type="checkbox"/>	Standard	Camber 3°		<input type="radio"/>
<input type="checkbox"/>		Camber 6°		<input type="radio"/>
Rear wheels type				
<input type="checkbox"/>	Standard	26" High Performance clincher wheels	With black spokes and alloy rims	<input type="radio"/>
<input type="checkbox"/>	XLT75	26" Top End Carbon rim / stainless spokes	Tubular only - Composite rims - Spokes: 24 front and 20 rear	£1,360
<input type="checkbox"/>	XLT78	26" Top End/Spinergy® Carbon rim	Clincher only	£1,288
<input type="checkbox"/>	XLT77	Time trial package: Zipp disc rear wheels / 26" Top End Carbon front wheel	ZIPP disc rear wheels are tubular only Front wheel with carbon rim	£2,555
<input type="checkbox"/>	XLT67	Corima® AERO 26" Carbon rim	Tubular only - With carbon/graphite composite rim and stainless ovalized spokes	£2,266

<input type="checkbox"/>	XLT62	Build for 700C wheels	Tubular only	£391
Tyres				
<input type="checkbox"/>	XLT51	Continental® Grand Prix tyres	Clincher only	£91
<input type="checkbox"/>	XLT52	Continental® Sprinter tyres	Tubular only	£84
SAFETY AND COMFORT OPTIONS				
	Standard	Seat and back cushions		○
	Standard	Drafting bumper, chain guard & safety flag		○
<input type="checkbox"/>	XLT71	Tool and tyre repair kit	Incl. 22 tools in 1, a unique 3 piece design plus tyre patch kit, CO2 tyre inflator, tyre changing tools and tool bag	£56
<input type="checkbox"/>	XLT16	Extra safety flag		£20
<input type="checkbox"/>	XLT63	Mirror, handle mount		£20
<input type="checkbox"/>	XLT13	Water bottle & cage		£24
<input type="checkbox"/>	XLT28	Backpack hydration system		£66
<input type="checkbox"/>	XLT141	Heart rate monitor		£66
<input type="checkbox"/>	XLT14	Cordless computer		£50
<input type="checkbox"/>	XLT15	Helmet		£50
	<input type="checkbox"/>	Small (6 5/8" - 7")		
	<input type="checkbox"/>	Medium (6 7/8" - 7 1/2")		
	<input type="checkbox"/>	Large (7 1/4 - 7 5/8")		
<input type="checkbox"/>	XLT9	Leg guard attachment		£50
<input type="checkbox"/>	XLT6	Front safety light		£17
<input type="checkbox"/>	XLT7	Rear flashing light		£17
<input type="checkbox"/>	XLT33	Alignment gauge		£50
<input type="checkbox"/>	XLT40	Cable kit		£25
<input type="checkbox"/>	XLT74	Extra steering dampener		£11
<input type="checkbox"/>	XLT42	Indoor training roller		£170
<input type="checkbox"/>	TER126	Click-strap large	Pad length approx. 12"	£58
<input type="checkbox"/>	TER125	Click-strap medium	Pad length approx. 10"	£54
<input type="checkbox"/>	TER124	Click-strap small	Pad length approx. 8"	£51
<input type="checkbox"/>	RACE4	Wheel bag		£180
<input type="checkbox"/>	RACE20	Handcycle travel bag		£300
<input type="checkbox"/>	RACE13	Top End bike jersey		£62
	<input type="checkbox"/>	L		
	<input type="checkbox"/>	XL		
<input type="checkbox"/>	RACE22	Top End t-shirt		£9
	<input type="checkbox"/>	L		
	<input type="checkbox"/>	XL		
UPHOLSTERY				
	Standard	Black breathable mesh upholstery with vertical neck support and 2" seat and back cushion standard. Seat restraint included.		○
FRAME COLOUR				
Contemporary colours				
<input type="checkbox"/>		Wet black		○
<input type="checkbox"/>		Black prism		○
<input type="checkbox"/>		White snow		○
<input type="checkbox"/>		Silver metallic		○
<input type="checkbox"/>		Stealth black		○
Way cool colours				
<input type="checkbox"/>		Royal blue		○
<input type="checkbox"/>		Midnight blue		○
<input type="checkbox"/>		Grasshopper green		○
Too hot colours				
<input type="checkbox"/>		Tangerine red		○
<input type="checkbox"/>		Raspberry red		○
<input type="checkbox"/>		Grape madness		○
<input type="checkbox"/>		Safety yellow		○
<input type="checkbox"/>		Sunkissed orange		○
<input type="checkbox"/>		Bubblegum pink		○
<input type="checkbox"/>	P1	Two colours option		£88
		Fork colour: _____		
		Frame colour: _____		

XLT81 Fork: Fixed height and length required measurements



Material information
Frame: Ovalized 7005 aluminium tubing with internal skeletal structure for increased stiffness. Wrap-around ovalized lower frame.
Fork: 7005 aluminium, super stiff with adjustable crank housing
Turning radius: 18 feet for 360° circle, 12 feet with one backing manoeuvre
Static tip angle: 39° for 250 lb/113Kg subject
Transfer seat height: 5" with 2" cushion (standard)
Backrest: adjustable fore and aft
Back angle: adjustable
Propulsion system: features a combination of SRAM and Top End components, 27 speeds, 48-39-24 chain rings and 11/34 cassette as standard. The 11/36 cassette is operated by the Hands-On Grip-Shifter or Trigger shifter unit which is mounted on the right (standard) hand pedal. It operates both the 11/34 cassette and the hands-on brake. A shifter/derailleur mounted on upper crank operates the top chain rings utilizing a system which minimizes cables and their maintenance. A parking brake/second brake is mounted on the fork.
Steering: pivot fork steering
Warranty: 3 year on frame and fork to original owner/limited parts warranty

Invacare Ltd: www.invacare.co.uk
 Pencoed Technology Park
 Pencoed
 Bridgend, CF35 5AQ T: 01656 776220 F: 01656 776222

All rights reserved. Trademarks are identified by the symbols ™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted.
 Shimano, Ultegra and Rapid Fire are registered trademarks of Shimano Inc.
 Spinerig is a registered trademark for Spinerig Corp., Corima is a registered trademark of Corima SA., Zipp is a registered trademark of Compositech, Inc.
 Cane Creek is a registered trademark of Cane Creek Cycling Components.
 SRM is a registered trademark of SRM Service Center Inc. USA
 SRAM is a registered trademark of SRAM Corp. Tufo is a registered trademark of Tufo of America Incorporated.
 Continental is a registered trademark of Continental Tire Corp. Rotor is a registered trademark of Rotor Components, LLC.
 ©2011 Invacare Corporation
 Form No. 09-013 Rev 07.11 - P2007-029

