

Invacare® Top End® Excelerator Handcycle-UK Price List and Order Form

For accuracy of ordering, please fax order to Invacare Customer Service on 01656 776220
or e-mail to: ordersuk@invacare.com

www.topendwheelchair.co.uk



Provider name: _____ Account number: _____
PO number: _____
Phone number: _____ Fax number: _____
Address: _____

Customer name: _____ Phone number: _____
Address: _____
Ship to: Provider Customer (Order will be shipped to provider unless otherwise noted)

DATE:

QUOTE


ORDER

- Please indicate your choice with a cross in this square
 Indicate that this option has no upcharge

Manual Sport Wheelchair
January 2012


Item	Description	Comments and/or non compatibilities information	Price
<input type="checkbox"/> XCL	Top End Excelerator		£1,899
FEATURES			
Standard	Carbon steel frame and fork		<input type="radio"/>
Standard	Adjustable seat and back height		<input type="radio"/>
Standard	Parking brake	Caliper brake with locking brake lever	<input type="radio"/>
Standard	9° camber	For stability	<input type="radio"/>
SEATING OPTIONS			
Seat width			
<input type="checkbox"/>	Seat width 14" (approx. 356 mm)		<input type="radio"/>
<input type="checkbox"/> XCL21	Custom: Seat width 15" (approx. 381 mm)		£58
<input type="checkbox"/>	Seat width 16" (approx. 406 mm)		<input type="radio"/>
<input type="checkbox"/> XCL21	Custom: Seat width 17" (approx. 432 mm)		£58
<input type="checkbox"/>	Seat width 18" (approx. 457 mm)		<input type="radio"/>
<input type="checkbox"/> XCL21	Custom: Seat width 19" (approx. 483 mm)		£58
<input type="checkbox"/>	Seat width 20" (approx. 508 mm)		<input type="radio"/>
FOOTREST			
Standard	Adjustable height/angle footrest system	Dual footrests with safety straps	<input type="radio"/>
PROPULSION			
<input type="checkbox"/> Standard	Shimano® 7 speed hub with reversing drum brake 175 mm length crankset with one chain ring	7 speeds, internal gearing, front wheel drive, 7 speed hub with internal hub brake and 7 speed twist shifter <u>Braking system:</u> Shimano internal drum/coaster brake activated by pedaling backwards. <u>Gear shifter:</u> Shimano Tourney revoshifter indexed twist grip with gear display. Check also quad options.	<input type="radio"/>
XCL2	Mountain drive	Only with std cranks and std propulsion system (7 speed internal hub). To gear "down" for climbing hills if you need it. Gives you a true doubling of the number of gears. Cannot change gears under load. Not available with V crankset option. Chain ring guard only.	£464
XLT139	Quad twist shifter option		£50
CRANK SET AND HAND PEDALS			
Standard	175 mm length crankset with one chain ring	Sugino 175 mm with chainring, (38 teeth)	<input type="radio"/>
Standard	Adjustable crank height		<input type="radio"/>
Standard	Chain tension idler	Included	<input type="radio"/>
Standard	Ovalized ergonomic handpedals	Shape and angle is ergonomically designed for max power and are a lightweight ovalized aluminium shape with frictionless bearings	<input type="radio"/>
<input type="checkbox"/> XLT8	Tri-pin quad handpedals	Tri-pin quad handles offer a 3 point contact (the same are used in driving a car) and also allow customer to get on/off of handles independently, as to change gears, etc.	£237
<input type="checkbox"/> XLT22	Quad gloves		
<input type="checkbox"/>	<input type="checkbox"/> Small <input type="checkbox"/> Medium <input type="checkbox"/> Large <input type="checkbox"/> X-Large	Quad gloves secure paraplegic's hand	£81
<input type="checkbox"/> XLT34	Top End V crankset with ergonomic handpedals Length: 180mm - Width: 22"	Not available with mountain drive option. Ergonomic aluminium ovalized handpedals are standard. Request to be made directly at Top End for quad hand pedal options.	£136

WHEELS				
<input type="checkbox"/>	Standard	Quick release axles	Easier for transport	<input type="radio"/>
<input type="checkbox"/>		Threaded axles	More performance. Screw into the receiver with allen key	<input type="radio"/>
Rear wheels type				
<input type="checkbox"/>	Standard	24" Cruiser wheels with cruiser tyres	Rear hubs are precision heat treated, silver anodized with quick release stainless axles or optional ½" threaded axles - 14 gauge stainless steel spokes.	<input type="radio"/>
Tyres				
<input type="checkbox"/>	Standard	Cruiser tyres	Cruiser wheels: 24" x 1.75 treaded tyre, tube 40 - 65 PSI	<input type="radio"/>
<input type="checkbox"/>	XCL22	Thorn resistant tyre protectors		£33
SAFETY AND COMFORT OPTIONS				
<input type="checkbox"/>	Standard	Seat cushion		<input type="radio"/>
<input type="checkbox"/>	Standard	Chain guard	Black plastic	<input type="radio"/>
<input type="checkbox"/>	Standard	Safety flag		<input type="radio"/>
<input type="checkbox"/>	Standard	Seat restraint: hook and loop		<input type="radio"/>
<input type="checkbox"/>	Standard	Reflectors		<input type="radio"/>
<input type="checkbox"/>	XLT13	Water bottle & cage		£24
<input type="checkbox"/>	XLT19	Tow bar for wheelchair		£170
<input type="checkbox"/>	XLT14	Cordless computer		£50
<input type="checkbox"/>	XLT6	Front safety light		£17
<input type="checkbox"/>	XLT7	Rear flashing light		£17
<input type="checkbox"/>	XLT15	Helmet		£50
<input type="checkbox"/>	<input type="checkbox"/> Small (6 5/8" - 7")			
<input type="checkbox"/>	<input type="checkbox"/> Medium (6 7/8" - 7 ½")			
<input type="checkbox"/>	<input type="checkbox"/> Large (7 ¼ - 7 5/8")			
<input type="checkbox"/>	XCL13	Seat restraint: auto style fastener		£33
<input type="checkbox"/>	TER126	Click-strap large	Pad length approx. 12" (305 mm)	£58
<input type="checkbox"/>	TER125	Click-strap medium	Pad length approx. 10" (254 mm)	£54
<input type="checkbox"/>	TER124	Click-strap small	Pad length approx. 8" (203 mm)	£51
<input type="checkbox"/>	XLT12	Mirror		£20
<input type="checkbox"/>	XLT17	Crutch holder and strap		£84
<input type="checkbox"/>	XLT33	Alignment gauge		£50
<input type="checkbox"/>	XLT42	Indoor training roller		£170
<input type="checkbox"/>	RACE22	Top End t-shirt		£9
<input type="checkbox"/>	<input type="checkbox"/> L			
<input type="checkbox"/>	<input type="checkbox"/> XL			
<input type="checkbox"/>	RACE13	Top End bike jersey		£62
<input type="checkbox"/>	<input type="checkbox"/> L			
<input type="checkbox"/>	<input type="checkbox"/> XL			
UPHOLSTERY				
<input type="checkbox"/>	Black Standard	Black nylon with hook and loop tension adjustable seat sling and seat back. Foam seat cushion with cover.		<input type="radio"/>
FRAME COLOUR				
Contemporary colours				
<input type="checkbox"/>		Wet black		<input type="radio"/>
<input type="checkbox"/>		Black prism		<input type="radio"/>
<input type="checkbox"/>		White snow		<input type="radio"/>
<input type="checkbox"/>		Silver metallic		<input type="radio"/>
<input type="checkbox"/>		Stealth black		<input type="radio"/>
Way cool colours				
<input type="checkbox"/>		Electric blue		<input type="radio"/>
<input type="checkbox"/>		Midnight blue		<input type="radio"/>
<input type="checkbox"/>		Grasshopper green		<input type="radio"/>
Too hot colours				
<input type="checkbox"/>		Tangerine red		<input type="radio"/>
<input type="checkbox"/>		Electric red		<input type="radio"/>
<input type="checkbox"/>		Grape madness		<input type="radio"/>
<input type="checkbox"/>		Safety yellow		<input type="radio"/>
<input type="checkbox"/>		Sunkissed orange		<input type="radio"/>
<input type="checkbox"/>		Bubblegum pink		<input type="radio"/>



Material information
Frame: 1 3/4" x .0065 ovalized carbon steel tubing
Frame design: single main tube, upright position
Wheelbase: overall length 65"
Weight: overall weight 50 lbs (approx. 22.7 kg)
Turning radius: 12 feet for 360° circle, 8 feet with one backing manoeuvre
Static tip angle: 25-30° for 250 lb / 113 kg subject
Transfer seat height: (with 2" cushion) 20" (approx. 508 mm)
Warranty: 3 years on frame and fork to original owner/limited parts warranty - Limited warranty on mountain drive (XCL2) option
Shipping: shipped completely assembled except for rear wheels, shipping weight is 80 lbs (approx. 36.3 kg)

Measurement guide for V crankset / handpedals



Handpedals are made of lightweight aluminium and are ovalized to precisely fit the hand. The bearing which allow the pedal to rotate while pedaling are located in the crank set and not in the pedal, which makes this pedal very comfortable. Includes one chainring. Mountain drive option not compatible with V Crankset option. See drawing for measurements. The width is measured horizontally centre to centre to the point at the middle of the hand pedal. The length is measured perpendicular from the centre of the crank to the centre of the hand pedal spindle.

Invacare Ltd: www.invacare.co.uk

Pencoed Technology Park, Pencoed, Bridgend, CF35 5AQ T: 01656 776222 F: 01656 776220

All rights reserved. Trademarks are identified by the symbols ™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted.

Shimano is a registered trademarks of Shimano Inc.

Delrin is registered trademark of E.I. duPont de Nemours and Company.

Specifications and prices are subject to change without notification.

©2011 Invacare Corporation - Form No. 96-266

