

# Invacare Top End Crossfire T6 Ultra Lightweight Active Performance Rigid Wheelchair - UK Price List and Order Form

For accuracy of ordering, please fax order to Invacare Customer Service on 01656 776220 or e-mail to: [ordersuk@invacare.com](mailto:ordersuk@invacare.com)

[www.topendwheelchair.co.uk](http://www.topendwheelchair.co.uk)



Provider name: \_\_\_\_\_ Account number: \_\_\_\_\_  
 PO number: \_\_\_\_\_  
 Phone number: \_\_\_\_\_ Fax number: \_\_\_\_\_  
 Address: \_\_\_\_\_

Customer name: \_\_\_\_\_ Phone number: \_\_\_\_\_  
 Address: \_\_\_\_\_

Ship to:  Provider  Customer (Order will be shipped to provider unless otherwise noted)

DATE:

QUOTE

ORDER

Manual Sport Wheelchair  
 January 2012

- Please indicate your choice with a cross in this square  
 Indicate that this option has no upcharge

Item	Description	Comments and/or non compatibilities information	Price
<input type="checkbox"/> CRF	Top End Crossfire T6	Aluminium wheelchair built of 6061 T6 Aluminium	£1,915
<b>FRONT FRAME MODIFICATIONS</b>			
<input type="checkbox"/> TFE	Tapered front end	N/C with front seat heights of 16" or 17"	<input type="radio"/>
<input type="checkbox"/> CRF2	Non tapered "V" front end	N/C with tapered seat width or tubular footrest	£134
<input type="checkbox"/> CRF1	Tapered seat width	N/C with ergo seat or widths of 12" or 13"	£102
<b>Taper difference rear to front (select one)</b>			
<input type="checkbox"/> CRFDRF1	Difference -1"	Tapers rear to front up to a maximum difference of 4". The front seat width measures the same as the outside of the footrest or straight so there is no tapered or V front end. Foot space is 2.25" less than front frame width.	<input type="radio"/>
<input type="checkbox"/> CRFDRF2	Difference -2"		<input type="radio"/>
<input type="checkbox"/> CRFDRF3	Difference -3"		<input type="radio"/>
<input type="checkbox"/> CRFDRF4	Difference -4"		<input type="radio"/>
<b>FORK OPTIONS</b>			
<input type="checkbox"/> CRF59	Frog legs® suspension forks	NC with 85° or 90° front frame angle or seat depths less than 16"	£201
<input type="checkbox"/> FLOSF	Frog legs® one sided forks	N/C with suspension forks / Black anodized colour only	£128
<b>SEATING OPTIONS</b>			
<b>Seat width</b>			
<input type="checkbox"/>	Seat width 12" (approx. 305 mm)		<input type="radio"/>
<input type="checkbox"/>	Seat width 13" (approx. 330 mm)		<input type="radio"/>
<input type="checkbox"/>	Seat width 14" (approx. 356 mm)		<input type="radio"/>
<input type="checkbox"/>	Seat width 15" (approx. 381 mm)		<input type="radio"/>
<input type="checkbox"/>	Seat width 16" (approx. 406 mm)		<input type="radio"/>
<input type="checkbox"/>	Seat width 17" (approx. 432 mm)		<input type="radio"/>
<input type="checkbox"/>	Seat width 18" (approx. 457 mm)		<input type="radio"/>
<input type="checkbox"/> CRF105	Seat width 19" (approx. 483 mm)		£67
<input type="checkbox"/> CRF105	Seat width 20" (approx. 508 mm)		£67
<b>Seat depth</b>			
<input type="checkbox"/>	Seat depth 14" (approx. 356 mm)		<input type="radio"/>
<input type="checkbox"/>	Seat depth 15" (approx. 381 mm)		<input type="radio"/>
<input type="checkbox"/>	Seat depth 16" (approx. 406 mm)		<input type="radio"/>
<input type="checkbox"/>	Seat depth 17" (approx. 432 mm)		<input type="radio"/>
<input type="checkbox"/>	Seat depth 18" (approx. 457 mm)		<input type="radio"/>
<input type="checkbox"/> CRF106	Seat depth 19" (approx. 483 mm)		£67
<input type="checkbox"/> CRF106	Seat depth 20" (approx. 508 mm)		£67
<b>Extra frame depth</b>			
<input type="checkbox"/>	0" (approx. 0 mm)		<input type="radio"/>
<input type="checkbox"/>	1" (approx. 25 mm)		<input type="radio"/>
<input type="checkbox"/>	1 ½" (approx. 38 mm)		<input type="radio"/>
<input type="checkbox"/>	2" (approx. 51 mm)		<input type="radio"/>
<b>SEAT HEIGHTS OPTIONS</b>			
<b>Seat height front</b>			
<input type="checkbox"/>	Seat height front 17" (approx. 432 mm)		<input type="radio"/>
<input type="checkbox"/>	Seat height front 17.5" (approx. 445 mm)		<input type="radio"/>
<input type="checkbox"/>	Seat height front 18" (approx. 457 mm)		<input type="radio"/>
<input type="checkbox"/>	Seat height front 18.5" (approx. 470 mm)		<input type="radio"/>
<input type="checkbox"/>	Seat height front 19" (approx. 483 mm)		<input type="radio"/>
<input type="checkbox"/>	Seat height front 19.5" (approx. 495 mm)		<input type="radio"/>
<input type="checkbox"/>	Seat height front 20" (approx. 508 mm)		<input type="radio"/>
<input type="checkbox"/>	Seat height front 20.5" (approx. 521 mm)		<input type="radio"/>
<input type="checkbox"/>	Seat height front 21" (approx. 533 mm)		<input type="radio"/>


Seat height rear				
<input type="checkbox"/>		Seat height rear 15" (approx. 381 mm)	Lowest rear seat height for 24" wheels with std tyres	<input type="radio"/>
<input type="checkbox"/>		Seat height rear 15.5" (approx. 394 mm)	Lowest rear seat height for 25" wheels with std tyres	<input type="radio"/>
<input type="checkbox"/>		Seat height rear 16" (approx. 406 mm)		<input type="radio"/>
<input type="checkbox"/>		Seat height rear 16.5" (approx. 419 mm)		<input type="radio"/>
<input type="checkbox"/>		Seat height rear 17" (approx. 432 mm)		<input type="radio"/>
<input type="checkbox"/>		Seat height rear 17.5" (approx. 445 mm)		<input type="radio"/>
<input type="checkbox"/>		Seat height rear 18" (approx. 457 mm)		<input type="radio"/>
<input type="checkbox"/>		Seat height rear 18.5" (approx. 470 mm)		<input type="radio"/>
<input type="checkbox"/>		Seat height rear 19" (approx. 483 mm)		<input type="radio"/>
<input type="checkbox"/>		Seat height rear 19.5" (approx. 495 mm)		<input type="radio"/>
<input type="checkbox"/>		Seat height rear 20" (approx. 508 mm)		<input type="radio"/>
<input type="checkbox"/>		Seat height rear 20.5" (approx. 521 mm)		<input type="radio"/>
<input type="checkbox"/>		Seat height rear 21" (approx. 533 mm)		<input type="radio"/>
ERGO SEAT OPTIONS				
<input type="checkbox"/>	CRF104	Ergo seat with tapered front end	N/C with tapered seat width. 15" is the lowest rear seat height for	£170
<input type="checkbox"/>	CRF110	Ergo seat with V style front end	24" wheels w/std tyres and 15.5" for 25" wheels w/std tyres	£240
FRONT FRAME ANGLE				
<input type="checkbox"/>	FFA70	70°		<input type="radio"/>
<input type="checkbox"/>	FFA75	75°		<input type="radio"/>
<input type="checkbox"/>	FFA80	80°		<input type="radio"/>
<input type="checkbox"/>	FFA85	85°		<input type="radio"/>
<input type="checkbox"/>	FFA90	90°		<input type="radio"/>
CENTRE OF GRAVITY (choose one)				
<input type="checkbox"/>	CRF111	Active position 1.25" - 5.5"	Measured from the front of the backpost to the centre of the rear axle	<input type="radio"/>
<input type="checkbox"/>	CRF113	Extended position -2.25" - 2.0"	N/C with tapered seat width	£120
CAMBER SYSTEM				
<input type="checkbox"/>	CRF630 (std)	0° camber tube with single camber		<input type="radio"/>
<input type="checkbox"/>	CRF633	3° camber tube with single camber		<input type="radio"/>
<input type="checkbox"/>	CRF636	6° camber tube with single camber		<input type="radio"/>
<input type="checkbox"/>	CRF639	9° camber tube with single camber		<input type="radio"/>
<input type="checkbox"/>	CRF640	0° extra camber tube		£100
<input type="checkbox"/>	CRF643	3° extra camber tube		£100
<input type="checkbox"/>	CRF646	6° extra camber tube		£100
<input type="checkbox"/>	CRF649	9° extra camber tube		£100
FOOT SPACE WIDTH				
	Standard	As per diagram	Exceptions: SW 12", 13", 14" have foot space of 8"	<input type="radio"/>
BACKREST				
Back type				
<input type="checkbox"/>	LKDWN (std)	Lock-down, folding ht/angle adj. Back	Manufactured of 1" aluminium tubing, only available with adjustable height (3/4" tube telescoping in 1" tube)	<input type="radio"/>
<input type="checkbox"/>	CRF1070	Fixed back with fixed angle welded at 0°, adj ht	Manufactured of 1" aluminium available in fixed or adj height (7/8" telescoping in 1" tube), Max height on fixed back is 16". Welded integral push handles are only available with fixed back height.	£133
<input type="checkbox"/>	CRF1073	Fixed back with fixed angle welded at 3°, adj ht		£133
<input type="checkbox"/>	CRF1076	Fixed back with fixed angle welded at 6°, adj ht		£133
<input type="checkbox"/>	CRF1079	Fixed back with fixed angle welded at 9°, adj ht		£133
<input type="checkbox"/>	Standard	Backrest crossbar opening (3.5" std)		<input type="radio"/>
<input type="checkbox"/>	TER196	Deep crossbar opening 4.5"		£30
<input type="checkbox"/>	TER196	Deep crossbar opening 5.5"		£30

<b>Back height</b>			
<input type="checkbox"/>	BH811	9"-11" (approx. 229 - 279 mm)	<input type="radio"/>
<input type="checkbox"/>	BH1014	10"-14" (approx. 254 - 356 mm)	<input type="radio"/>
<input type="checkbox"/>	BH1216	12"-16" (approx. 305 - 406 mm)	<input type="radio"/>
<input type="checkbox"/>	BH1418	14"-18" (approx. 356 - 457 mm)	<input type="radio"/>
<input type="checkbox"/>	BH1620	16"-20" (approx. 406 - 508 mm)	<input type="radio"/>
<b>Back upholstery and accessories</b>			
<input type="checkbox"/>	CRFBAT	Ballistic nylon ten adj with straps and anchors	<input type="radio"/>
<input type="checkbox"/>	MESH	Breathable mesh tension adj with anchors	<input type="radio"/>
<input type="checkbox"/>	TER195	Back foam upgrade 1" extra	£126
<input type="checkbox"/>	TER197	Back centre width upgrade (equal to seat width)	£126
<input type="checkbox"/>	CRF16	Fold over back upholstery	N/C with integral (CRF97) push handles and camouflage uphol
<input type="checkbox"/>	CRF30	Back pack	£32
<input type="checkbox"/>	CRF82	Push handles height adjustable bolt-on	Bolted to the back post and have 10" of adjustment up or down. 16" tall. N/C with swingaway arms. £133
<input type="checkbox"/>	CRF13	Push handles non-adjustable bolt-on	£50
<input type="checkbox"/>	CRF97	Bolt-through integral push handles	£68
<b>WHEELS</b>			
<b>Front wheels / Forks</b>			
<input type="checkbox"/>	CRFF3	3" Micro	<input type="radio"/>
<input type="checkbox"/>	3SOFR	3" Aluminium soft roll	£60
<input type="checkbox"/>	CRF654	4" Semi-pneumatic	<input type="radio"/>
<input type="checkbox"/>	CRF655	5" Semi-pneumatic	N/C with 90° front frame angle <input type="radio"/>
<input type="checkbox"/>	CRFA654	4" Aluminium soft roll	Alu castors are 1.4" with soft roll tyres. For suspension, see p.1 £60
<input type="checkbox"/>	CRFA655	5" Aluminium soft roll	N/C with 90° front frame angle, 1.4" wide with soft roll tyres £60
<input type="checkbox"/>	CRF884	4" Lighted front wheels	£47
<b>Rear wheel options</b>			
<input type="checkbox"/>	CRF22SP	22" Spoke	<input type="radio"/>
<input type="checkbox"/>	CRF24HPS	24" Sun® High Performance spoke	<input type="radio"/>
<input type="checkbox"/>	CRF25HPS	25" Sun® High Performance spoke	<input type="radio"/>
<input type="checkbox"/>	CRF26HPS	26" Sun® High Performance spoke	£133
<input type="checkbox"/>	FN1624	24" Sun® Sun Fusion 16 spoke wheel	N/C with CRF73 £223
<input type="checkbox"/>	FN1625	25" Sun® Sun Fusion 16 spoke wheel	N/C with CRFPN and CRF73 £223
<input type="checkbox"/>	FN1626	26" Sun® Sun Fusion 16 spoke wheel	N/C with CRF73, CRF60, CRF7G, CRF7BK, CRF101R and CRF101Y £223
<input type="checkbox"/>	CRF6722	22" Spinergy SPOX wheels	N/C with CRF123, CRF60, CRF7G, CRF7BK, CRF101R and CRF101Y £328
<input type="checkbox"/>	CRF6724	24" Spinergy SPOX wheels	N/C with CRF73, CRF6725 £328
<input type="checkbox"/>	SPLX24	24" Spinergy LX wheels (12 spoke)	N/C with CRF73, CRF6725 £346
<input type="checkbox"/>	CRF6725	25" Spinergy SPOX wheels	N/C with CRFPN and CRF73 £328
<input type="checkbox"/>	SPLX25	25" Spinergy LX wheels (12 spoke)	N/C with CRFPN and CRF73 £346
<input type="checkbox"/>	CRF6726	26" Spinergy SPOX wheels	N/C with CRF73, CRF60, CRF7G, CRF7BK, CRF101R and CRF101Y £328
<b>Spoke colour (only for Spinergy® wheels)</b>			
<input type="checkbox"/>	SPCBK - Std	Black	If colour is not chosen, spokes will be black <input type="radio"/>
<input type="checkbox"/>	SPCY	Yellow	<input type="radio"/>
<input type="checkbox"/>	SPCBW	White	<input type="radio"/>
<input type="checkbox"/>	SPCR	Red	<input type="radio"/>
<input type="checkbox"/>	SPCBL	Blue	<input type="radio"/>
<input type="checkbox"/>	CRF27	Quad release axles	£84
<input type="checkbox"/>	CRF76	Titanium axles	£102
<b>Tyre options</b>			
<input type="checkbox"/>	CRFPR	Primo high pressure tyres: grey	<input type="radio"/>
<input type="checkbox"/>	CRFPN	Pneumatic - 1 3/8" treaded: grey	<input type="radio"/>
<input type="checkbox"/>	CRF123	Schwalbe 1" flat resistant: black	£60
<input type="checkbox"/>	CRF73	Pneu w/FF - 22 or 24x 1 3/8" tyres	£84
<input type="checkbox"/>	CRF60	Knobby tyres: grey	£50
<input type="checkbox"/>	CRF7G	Solid no flat tyres: grey	£50
<input type="checkbox"/>	CRF7BK	Solid no flat tyres: black	£50
<input type="checkbox"/>	CRF101R	Kenda® high pressure tyres: red	£50
<input type="checkbox"/>	CRF101Y	Kenda® high pressure tyres: yellow	£50
<b>ANODIZING ACCENTS</b>			
<input type="checkbox"/>	ANOBLK	Bold black (standard)	Includes standard forks and back angle bracket <input type="radio"/>
<input type="checkbox"/>	ANOSLV	Classic silver	Includes standard forks and back angle bracket <input type="radio"/>
<input type="checkbox"/>	TED140	Shadow blue	Includes standard forks and back angle bracket £50
<input type="checkbox"/>	TED198	Victory red	Includes standard forks and back angle bracket £50
<b>WHEELS: EXTRA SET FOR OFF ROAD (includes wheels, tyres and aluminium pushrims)</b>			
<input type="checkbox"/>	CRF91	24" Knobby tyres (540)	24" Knobby match up with 25" daily use wheels, with Primo tyres £374
<input type="checkbox"/>	CRF92	25" Knobby tyres (559)	25" Knobby match up with 26" daily use wheels, with Primo tyres £374
<input type="checkbox"/>	CRF93	Extra set of axles	£36

<b>HANDRIMS</b>			
<b>Mounting position</b>			
<input type="checkbox"/>	OUT	Standard mount	<input type="radio"/>
<input type="checkbox"/>	IN	Close mount	<input type="radio"/>
<b>Handrim options</b>			
<input type="checkbox"/>	CRF77 (std)	Aluminium anodized handrims	<input type="radio"/>
<input type="checkbox"/>	CRF75	Titanium handrims	N/C with 22" wheels £260
<input type="checkbox"/>	CRF8	Coated handrims	£71
<input type="checkbox"/>	NFSDLT	Natural Fit® LT standard grip handrims	N/C with 22" wheels £205
<input type="checkbox"/>	NFSPLT	Natural Fit® LT super grip handrims	N/C with 22" wheels £233
<input type="checkbox"/>	CRF99	Natural Fit® standard grip handrims	N/C with 22" wheels £187
<input type="checkbox"/>	CRF98	Natural Fit® super grip handrims	N/C with 22" wheels £217
<b>FOOTPLATE OPTIONS</b>			
<b>Type of footplate</b>			
<input type="checkbox"/>	CRFTFR	1 piece tubular footrest	Without cover <input type="radio"/>
<input type="checkbox"/>	CRF1FP	Adjustable angle one piece footplate	Standard for frames with "V" style front <input type="radio"/>
<input type="checkbox"/>	CRF191	Adjustable angle footrest with 2" risers	£64
<b>Footrest / Legrest accessories</b>			
<input type="checkbox"/>	CRF10	Plastic footrest cover	For tubular footrest only £50
<input type="checkbox"/>	CRF89	Carbon fibre footrest cover	For tubular footrest only £84
<input type="checkbox"/>	CRF28	Frame protector set (2 front, 1 rear)	£39
<input type="checkbox"/>	CRF103	Luggage carriers	£68
<input type="checkbox"/>	CRF133	Calf strap (standard)	<input type="radio"/>
<b>SEATING OPTIONS</b>			
<b>Seat accessories</b>			
<input type="checkbox"/>	CRF29	Seat pouch	£23
<input type="checkbox"/>	CRF72	Positioning strap - Hook & Loop	£26
<input type="checkbox"/>	CRF74	Positioning strap - Auto style	£33
<b>Cushions</b>			
<input type="checkbox"/>	CRF14	2" soft cushion	Filled with 2" soft foam which offer minimum protection £50
<b>Seat upholstery</b>			
<input type="checkbox"/>	CRFATS	Adj tension seat with straps & anchor loop (standard)	<input type="radio"/>
<input type="checkbox"/>	CRFT112	Screw-on seat upholstery	N/C with tapered seat width £52
<b>ARM OPTIONS</b>			
<b>Arm type</b>			
<input type="checkbox"/>	CRF78	Swingaway arms	Adjustable in 1" increments from 10.5" to 12.5" (top of the seat rail to top of armrest). N/C with adj height push handles. £141
<b>Arm accessories</b>			
<input type="checkbox"/>	CRF121	Fold-in style black aluminium side guards	Fold-in style flip back for easier transfers and are adjustable in height. Compatible with both high (push or pull to lock) or £84
<input type="checkbox"/>	CRFSG	Fold-in style carbon fibre side guards	undermount (scissor) wheel locks. NC with adj height push handles if back height is 10-14" or less. £203
<input type="checkbox"/>	CRF79	Removable side guards / fender	£120
<input type="checkbox"/>	CRF19	Removable side guards without hole	£84
<input type="checkbox"/>	CRF18	Removable side guards with hole (quad)	£84
<input type="checkbox"/>	CRF20E	Fixed aluminium side guards even w/wheels	Fixed back only £122
<input type="checkbox"/>	CRF201	Fixed aluminium side guards 1" above wheels	Fixed back only £122
<input type="checkbox"/>	CRF83E	Fixed carbon fibre side guards even w/wheels	Fixed back only, N/C with tapered seat width £170
<input type="checkbox"/>	CRF831	Fixed carbon fibre side guards 1" above wheels	Fixed back only, N/C with tapered seat width £170
<b>WHEEL LOCKS (choose one)</b>			
<input type="checkbox"/>	CRFPUL	Pull to lock wheel locks	<input type="radio"/>
<input type="checkbox"/>	CRFPSH	Push to lock wheel locks	<input type="radio"/>
<input type="checkbox"/>	CRF109	Hideaway undermount wheel locks	N/C with 14" depth £45
<input type="checkbox"/>	CRF22	Hill holders	£102
<input type="checkbox"/>	CRFOWL	Omit wheel locks	<input type="radio"/>
<b>MISCELLANEOUS ACCESSORIES</b>			
<input type="checkbox"/>	FRMPRBK	Chair armour black neoprene frame protectors	£39
<input type="checkbox"/>	FRMPRBK	Chair armour blue neoprene frame protectors	£39
<input type="checkbox"/>	FRMPRR	Chair armour red neoprene frame protectors	£39
<input type="checkbox"/>	CRF21C	Spoke protectors: clear	£52
<input type="checkbox"/>	CRF23	Rear anti tippers	£68
<input type="checkbox"/>	CRF94	Alignment gauge	£50

UPHOLSTERY COLOUR			
<input type="checkbox"/>	Standard	Black	○
<input type="checkbox"/>	UPHI	Camouflage	NC with screw-on seat, mesh or fold over back upholstery. Ballistic nylon tension adjustable with anchors and straps for the back upholstery. If cushion is ordered it will also be camouflage. £67
FRAME COLOUR			
Contemporary colours			
<input type="checkbox"/>		Wet black	○
<input type="checkbox"/>		Black prism	○
<input type="checkbox"/>		White snow	○
<input type="checkbox"/>		Silver metallic	○
<input type="checkbox"/>		Stealth black	○
Way cool colours			
<input type="checkbox"/>		Electric blue	○
<input type="checkbox"/>		Midnight blue	○
<input type="checkbox"/>		Grasshopper green	○
Too hot colours			
<input type="checkbox"/>		Tangerine red	○
<input type="checkbox"/>		Electric red	○
<input type="checkbox"/>		Grape madness	○
<input type="checkbox"/>		Safety yellow	○
<input type="checkbox"/>		Sunkissed orange	○
<input type="checkbox"/>		Bubblegum pink	○

**ADDITIONAL INFORMATION & MEASURING GUIDE**



Knee-to-Heel Range			
Front Seat-to-Floor			
17"	12.25" - 16.25"	20"	15.25" - 19.25"
18"	13.25" - 17.25"	21"	16.25" - 20.25"
19"	14.25" - 18.25"		

Wheel/Tire Chart – Applies to Spoke, Sun & Spinerge F/Wheels				
Tire Choices	22"	24"	25"	26"
Primo: 1" high pressure grey clincher	X	X	X	X
Kenda: 1" high pressure clincher		X	X	X
Schwalbe 1" Black flat resistant		X	X	X
1-3/8" grey treaded	X	X		X
1-3/8" w/ flat free	X	X		
Solid 1" no-flat		X	X	
Knobby: 2" wide grey		X	X	

Overall Width by Camber Angle				
Wheel Type	0°	3°	6°	9°
Spoke	6.75"	8"	10"	12.25"

Add the amount above to the seat width to obtain the overall width at its widest point.



Width & Depth Matrix									
	W12	W13	W14	W15	W16	W17	W18	W19	W20
D12	18.63	18.73	18.83	18.93	19.03	19.13	19.23	19.33	19.43
D13	18.73	18.83	18.93	19.03	19.13	19.23	19.33	19.43	19.53
D14	18.83	18.93	19.03	19.13	19.23	19.33	19.43	19.53	19.63
D15	18.93	19.03	19.13	19.23	19.33	19.43	19.53	19.63	19.73
D16	19.03	19.13	19.23	19.33	19.43	19.53	19.63	19.73	19.83
D17	19.13	19.23	19.33	19.43	19.53	19.63	19.73	19.83	19.93
D18	19.23	19.33	19.43	19.53	19.63	19.73	19.83	19.93	20.03
D19	19.33	19.43	19.53	19.63	19.73	19.83	19.93	20.03	20.13
D20	19.43	19.53	19.63	19.73	19.83	19.93	20.03	20.13	20.23

**Section 1**

**Standard Chair Includes:**

- Tubular footrest without plate
- 4" Semi-pneumatic casters/ forks
- 3° Camber
- 24" High performance wheels
- 24" High pressure tires
- Medium height folding back
- Adjustable tension seat sling
- Adjustable tension back upholstery

**Top End Crossfire Component Weights – Section 2**

<b>Seating</b>	
Omit adjustable tension seat upholstery	-.85
No cushion	0
2" cushion	+1.10
3" cushion	+1.65
Screw-on seat upholstery	+0.5
<b>Back Type</b>	
Folding	0
Fixed	-.80
<b>Push Handles</b>	
Adjustable height bolt-on	+1.60
Fixed height bolt-on	+1.70
Integral: welded	+0.75
Integral: bolt-through	+0.75
<b>Back Upholstery</b>	
Adjustable tension w/ straps	0
Fold over	+0.09
Omit back upholstery	-.97
<b>Armrests</b>	
Swingaway padded tubular	+2.64
<b>Footrest</b>	
Tubular	0
Adjustable angle one piece footplate	+0.60
Plastic footrest cover	+0.35
Carbon fiber footrest cover	-.30
<b>Front Casters (pair)</b>	
3" Micro	-.13
3" Softroll	+0.19
3" Micro lighted	-.13
4" Semi pneumatic	0
5" Semi pneumatic	+0.31
4" Aluminum softroll	+0.44
5" Aluminum softroll	+0.92
4" Lighted	+0.22
5" Lighted	+0.53
Frog leg front suspension forks	+0.64
Frog Leg one-sided forks	-.05
<b>Wheel Locks</b>	
Push or pull-to-lock	+1.00
Undermount scissor	+1.25
Hill holders/grade aids	+0.51
<b>Sideguards</b>	
Fold-in Style carbon fiber	+0.95

Fold-in Style black aluminum	+1.45
Removable	+1.45
Removable side guards with fender	+1.94
Fixed aluminum	+1.45
Fixed carbon fiber	+1.00
<b>Rear Tires (pair)</b>	
High pressure, Primo or Kenda	0
Schwalbe	+0.75
Pneumatic 1 3/8"	+1.40
Pneumatic Flat Free	+1.45
Knobby	+3.10
Solid No Flat	+2.15
<b>Rear Wheels (pair)</b>	
22" spoke	-.60
24" Sun high performance spoke	0
25" Sun high performance spoke	+0.30
26" Sun high performance spoke	+0.65
24" Spinerge Spox/Sun Fusion	-.80
24" Spinerge LX	-1.20
25" Spinerge Spox/Sun Fusion	-.60
25" Spinerge LX	-1.00
26" Spinerge Spox/Sun Fusion	-.25
<b>Handrim</b>	
Aluminum	0
Titanium	+0.20
Plastic coated	+1.40
Natural Fit standard or super grip	+1.70
Natural Fit LT	+1.25
<b>Axles</b>	
Standard quick release	0
Titanium quick release	+0.17
Quad quick release	+0.22
<b>Accessories</b>	
Chair armor	+0.10
Luggage carriers	+0.75
Calf strap	+0.10
Seat pouch	+0.35
Positioning Strap – hook and loop	+0.10
Positioning Strap – auto style	+0.50
Spoke guards – black	+1.40
Spoke guards – clear	+1.60
Anti-tip bars	+1.63
Crutch holder	+1.45

• Total Weight Section 1 \_\_\_\_\_ Total Weight \_\_\_\_\_  
 • Total Weight Section 2 \_\_\_\_\_

Invacare Ltd: [www.invacare.co.uk](http://www.invacare.co.uk)

Invacare Ltd Pencoed Technology Park, Pencoed, Bridgend, CF35 5AQ T: 01656 776222 F: 01656 776220

All rights reserved. Trademarks are identified by the symbols ™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted.

Sun® is a registered trademark of Sun Components Inc. Natural Fit is a registered trademark of Three Rivers Holding, LLC. Frog Legs is a registered trademark of Frog Legs, Inc.

Spinerge is a registered trademark for Spinerge Corp. Schwalbe is a registered trademark of RALF BOHLE GmbH in Reichshof, Germany.

Kenda is a registered trademark of Kenda Rubber Industrial Co. Ltd. Specifications and prices are subject to change without notification.

©2011 Invacare Corporation - Form No. 06-001

